



# Product and solutions catalog

Sleep and home respiratory care

**PHILIPS**

**RESPIRONICS**

# Table of contents

<b>Sleep therapy</b> .....	3
<b>Home ventilaton</b> .....	9
<b>Airway clearance</b> .....	16
<b>Respiratory drug delivery</b> .....	19
<b>Oxygen</b> .....	36
<b>Masks</b> .....	51
<b>Sleep diagnostics</b> .....	69
<b>Pro-tech sensors</b> .....	90
<b>Value-added programs</b> .....	107
<b>Medical education</b> .....	110



We're working with you to establish healthier business by streamlining patient management, increasing efficiency, and maximizing value.



# Sleep therapy

Philips Respironics provides a family of System One solutions designed to deliver exceptional therapy, enhance patient comfort, and provide essential compliance tools.

Not all products are available in all countries. Check with Philips Respironics customer service or your local distributor.



# System One: a quality choice

With your input, in addition to the enhanced features and appearance of System One, we've included even more **quality and durability**, while providing you the same exceptional product and sales support.

## System One quality comes through

- 25% quieter than its closest competitor
- Designed to eliminate water ingress
- Durable, easy-to-connect, top-mount tubing connections for an exceptional customer experience
- Commitment to quality excellence that earned Philips Respironics the prestigious 2011 Manufacturing Plant of the Year by *Assembly Magazine*
- Assembled in the USA

**Learn more.**







The SleepMapper self-management solution is compatible with most System One devices, including bi-level systems. To find out more, go to: [www.sleepmapper.com/compatible](http://www.sleepmapper.com/compatible)

## Your System One CPAP options

### REMstar Auto<sup>†</sup>

- 4 modes (CPAP, Auto, CPAP-Check, Auto-Trial)
- A-Flex for comfort
- Opti-Start technology
- Advanced event detection
- Heated Tube option

### REMstar Pro<sup>†</sup>

- 3 modes (CPAP, CPAP-Check, Auto-Trial)
- C-Flex+ for comfort
- Advanced event detection
- Heated Tube option

### REMstar Plus<sup>†</sup>

- 1 mode (CPAP)
- C-Flex for comfort
- Patterns of use data reporting
- Heated Tube option

### REMstar SE<sup>†</sup>

- Flex for comfort
- Enhanced VIC reporting
- On-board memory
- Heated Tube option

<sup>†</sup> Modem and SD card capable; optional integrated oximetry available.

Ordering information	US part number	Canadian part number	HCPCS
<b>REMstar Auto</b>			
REMstar Auto A-Flex	DS560S	CA561S	E0601
REMstar Auto A-Flex with humidifier	DS560HS	CA561HS	E0601, E0562
REMstar Auto A-Flex with Heated Tube humidifier	DS560TS	CA561TS	E0601, E0562, A4604
<b>REMstar Pro</b>			
REMstar Pro C-Flex+	DS460S	CA461S	E0601
REMstar Pro C-Flex+ with humidifier	DS460HS	CA461HS	E0601, E0562
REMstar Pro C-Flex+ with Heated Tube humidifier	DS460TS	CA461TS	E0601, E0562, A4604
<b>REMstar Plus</b>			
REMstar Plus C-Flex	DS260S	CA261S	E0601
REMstar Plus C-Flex with humidifier	DS260HS	CA261HS	E0601, E0562
REMstar Plus C-Flex with Heated Tube humidifier	DS260TS	CA261TS	E0601, E0562, A4604
<b>REMstar SE</b>			
REMstar SE	DS220S	CA221S	E0601
REMstar SE with humidifier	DS220HS	CA221HS	E0601, E0562
REMstar SE with Heated Tube humidifier	DS220TS	CA221TS	E0601, E0562, A4604

### Product warranty

All System One devices and humidifiers come standard with a two-year warranty. US only.

# Therapy for your noncompliant, nontolerant CPAP patients

## Your System One bi-level options

Philips Respironics revolutionized sleep therapy by introducing bi-level positive airway pressure technology to treat obstructive sleep apnea. Emulating natural breathing with pressure relief at inhalation and exhalation, this proven technology does more than provide highly effective therapy – it delivers a level of comfort that can mean the difference between compliance and non-compliance for patients having difficulty adapting to traditional CPAP therapy.

### BiPAP Auto

- 3 modes (CPAP, fixed bi-level, auto bi-level)
- Adjustable P<sub>S</sub>min
- Bi-Flex for comfort
- Advanced event detection
- Heated Tube option
- Modem and SD card capable
- Optional integrated oximetry available

### BiPAP Pro

- 2 modes (CPAP, fixed bi-level)
- Bi-Flex for comfort
- Advanced event detection
- Heated Tube option
- Modem and SD card capable
- Optional integrated oximetry available



The SleepMapper self-management solution is compatible with most System One devices, including bi-level systems. To find out more, go to: [www.sleepmapper.com/compatible](http://www.sleepmapper.com/compatible)

Ordering information	US part number	Canadian part number	HCPCS
<b>BiPAP Auto</b>			
BiPAP Auto Bi-Flex	DS760S	CA761S	E0470
BiPAP Auto Bi-Flex with humidifier	DS760HS	CA761HS	E0470, E0562
BiPAP Auto Bi-Flex with Heated Tube humidifier	DS760TS	CA761TS	E0470, E0562, A4604
<b>BiPAP Pro</b>			
BiPAP Pro Bi-Flex	DS660S	CA661S	E0470
BiPAP Pro Bi-Flex with humidifier	DS660HS	CA661HS	E0470, E0562
BiPAP Pro Bi-Flex with Heated Tube humidifier	DS660TS	CA661TS	E0470, E0562, A4604

### Product warranty

All System One devices and humidifiers come standard with a two-year warranty. US only.

# Therapy for your most complicated patients

## Your System One advanced therapy options

### BiPAP autoSV Advanced

Targets: complex, central, mixed sleep apneas, and Cheyne-Stokes respiration patients

- Clinically proven servo ventilation algorithm
- Auto EPAP for a patient's changing needs
- Real time access to patient data via modems
- Advanced event detection
- Bi-Flex for comfort
- Encore data management and reporting
- New patient waveform reports
- Optional integrated oximetry available



Ordering information	US part number	Canadian part number	HCPCS
BiPAP autoSV Advanced	DS960S	CA961S	E0471
BiPAP autoSV Advanced with humidifier	DS960HS	CA961HS	E0471, E0562
BiPAP autoSV Advanced with Heated Tube humidifier	DS960TS	CA961TS	E0471, E0562, A4604

### BiPAP AVAPS

Targets: COPD, obesity hypoventilation, restrictive disorders, and neuromuscular disorders

- AVAPS feature automatically adjusts pressure support to maintain a target tidal volume
- Automatically adapts to disease progression and changing patient needs
- Bi-Flex for comfort
- Improves ventilation efficacy<sup>1</sup>
- Encore data management and reporting
- Digital Auto-Trak provides advanced leak compensation and automatic triggering and cycling
- Modem and SD card capable
- Optional integrated oximetry available



Ordering information	US part number	Canadian part number	HCPCS
BiPAP S/T, C Series	DS1060S	NA	E0471
BiPAP S/T, C Series, Heated Tube core package	DS1060TS	NA	E0471, E0562, A4604
BiPAP AVAPS, C Series	DS1160S	CA1161X	E0471
BiPAP AVAPS, C Series, Heated Tube core package	DS1160TS	CA1161XTS	E0471, E0562, A4604

### Product warranty

All System One devices and humidifiers come standard with a two-year warranty. US only.

<sup>1</sup> Storre, JH, et al. Average Volume-Assured Pressure Support in Obesity Hypoventilation: A Randomized Crossover Trial. Chest 2006;130:815-821



# Sleep therapy accessories

Sleep therapy accessories		
Humidification	Part number	HCPCS
System One 60 Series humidifier	DS6H, IN6H (Canada)	E0562
System One 60 Series Heated Tube humidifier	DS6T, IN6T (Canada) (Includes 80 watt power supply and Heated Tube)	E0562, A4604
Water chamber (includes both 50 and 60 Series lids)	1063785	A7046
Heated Tube humidifier upgrade kit	6T-CONVKIT	A4604
Power supplies		
60 watt power supply (for H and HS models)	1091398	
80 watt power supply (for TS model only)	1091399	
Power cord	1005894	
Tubing		
System One Heated Tube	SYS1HT15	A4604
System One Performance Tube (15mm)	1074400	A7037
Performance Tubing (22mm)	1032907	A7037
Filters		
Pollen (foam) filter, reusable (1 per pack)	1035443	A7039
Pollen (foam) filter, reusable (2 per pack)	1029330	A7039
Ultra-fine filter, disposable (1 per pack)	1063096	A7038
Ultra-fine filter, disposable (2 per pack)	1035442	A7038
Ultra-fine filter, disposable (6 per pack)	1029331	A7038
Accessories		
Device side panel	1099592	
Carrying case	1063857	
Data management accessories		
Data cards and accessories		
SD card (10 pack)	1063859	
SD card mailer (10 pack)	1065146	
SD card cover (attached to device)	1063858	
SD card reader	1047300	
Modems		
Wired flow modem	100200W	A9279
Wired flow modem monthly service	100201W	
Wireless flow modem	100200C, 100210C (Canada)	
Wireless flow modem monthly service	100201C	
Broadband modem	100200E	A9279
Modem oximetry accessories		
Oximetry accessory	100205E	
Masimo HPLP oximetry module	1075383	
Masimo oximetry sensor	Various part numbers	
Travel accessories		
Shielded DC cord	1097568	
DC battery adapter cable	532209	
Shielded 12VDC power cord system (includes shielded DC cord and battery adapter cable)	1099011	
Universal portable battery pack	1028869	



# Home ventilation

Lightweight, versatile, and easy to use, Philips Respironics ventilation technologies are designed for use in the home, hospital and alternative care sites and provide invasive and noninvasive ventilator support for a wide range of adult and pediatric patients.

Not all products are available in all countries. Check with Philips Respironics customer service or your local distributor.

## BiPAP AVAPS

Part number		HCPGS
DS1160S	BiPAP AVAPS	E0471
DS1160HS	BiPAP AVAPS with humidifier	E0471, E0562
DS1160TS	BiPAP AVAPS, Heated Tube humidifier	E0471, E0562, A4604
CA1161X (Canada)	BiPAP AVAPS	
CA1161XTS (Canada)	BiPAP AVAPS Heated Tube core package	

To learn more visit [www.philips.us/bipapavaps](http://www.philips.us/bipapavaps)



### BiPAP AVAPS

Aligned with System One sleep therapy technology, the BiPAP AVAPS noninvasive ventilator provides the comfort of pressure ventilation and the consistent efficacy of assured tidal volume. Adapting automatically to changing patient needs, BiPAP AVAPS incorporates advanced technologies that help simplify patient care.



Part number		HCPCS
DS1060S	BiPAP S/T	E0471
DS1060HS	BiPAP S/T with humidifier	E0471, E0562
DS1060TS	BiPAP S/T with Heated Tube humidifier	E0471, E0562, A4604

Not available in Canada

Accessories		
Humidification		
DS6H	System One 60 Series humidifier	E0562
IN6H (Canada)		
DS6T	System One 60 Series Heated Tube humidifier (includes 80 watt power supply and Heated Tube)	E0562, A4604
IN6T (Canada)		
1063785	Water chamber (includes both 50 and 60 Series lids) replacement only	A7046
6T-CONVKIT	Heated Tube humidifier upgrade kit	A4604 – Heated Tubing
Tubing		
SYS1HT15	System One Heated Tube	A4604
1074400	System One performance tubing (15mm)	A7037
1032907	System One performance tubing (22mm)	A7037
Filters		
1035443	Pollen (foam) filter, reusable (1 pk)	A7039
1029330	Pollen (foam) filter, reusable (2 pk)	A7039
1063096	Ultra-fine filter, disposable (1 pk)	A7038
1035442	Ultra-fine filter, disposable (2 pk)	A7038
1029331	Ultra-fine filter, disposable (6 pk)	A7038

Accessories	
1099592	Device side panel
1063857	Carrying case



**BiPAP S/T**

The BiPAP S/T noninvasive ventilator offers a straightforward user interface and integrated heated humidification for easy setup and management. Digital Auto-Trak automatically adjusts variable trigger cycle thresholds enabling optimum patient/ventilator synchrony while System One Heated Tube humidification allows humidity levels to be tailored to meet patient needs.

To learn more visit [www.philips.us/bipapst](http://www.philips.us/bipapst)

# Market-leading versatility

Technologies that automatically adjust to changing patient needs



Lightweight, versatile and easy to use, Trilogy ventilators are designed for use in the home, hospital and alternative care sites and provide invasive and noninvasive ventilator support for a wide range of pediatric and adult patients. They combine Auto-Trak and AVAPS technologies with advanced flow triggering and leak compensation features that allow for a more sensitive delivery of therapy. Trilogy also incorporates mouthpiece ventilation (MPV). A 'kiss' trigger detects when a patient engages and disengages from the mouthpiece to deliver on-demand ventilation. It combines with an optional MPV support system to enhance patient ease of use.

#### Key features

- Mouthpiece ventilation (MPV)
- Pressure support up to 40 cm H<sub>2</sub>O
- AVAP-AE
- Adjustable AVAPS
- Sensitive Auto-Trak
- On-screen waveforms

To learn more visit [www.philips.com/trilogy100](http://www.philips.com/trilogy100)

## Trilogy100 and Trilogy200

Part number		HPCPS
1054260	Trilogy100 ventilator	E0463 or
CA1054096	Trilogy100 ventilator, Canada	E0464
1040005	Trilogy200 ventilator	E0463 or
CA1032800	Trilogy200 ventilator, Canada	E0464
Power cords		
1044760	Power cord, C7-end, RA, FemEnd, 8ft, USA/Can	
1069219	Power cord, C7-end, RA, FedEnd, 14 ft, USA/Can	
Cables		
1045290	Nurse call adaptor cable	
1069979	1/4" to 2.5mm NC nurse call cable	
1069980	1/4" to 2.5mm NO nurse call cable	
1045291	Remote alarm adapter cable	
1046972	Trilogy RS232 communication cable	
1080247	RJ9 to 1/4" nurse call cable alarm closed	
1080248	RJ9 to 1/4" nurse call cable alarm open	
1080249	Nurse call RJ9 cable open ends	
1080250	RJ9 to 2.5mm 3" female adapter cable	
1047295	External battery cable	A4612
1067424	Automotive DC power adapter	A4612
Porting blocks		
1054670	Trilogy100, exhalation port block, active	
1040372	Trilogy100, exhalation port block, passive	
1040370	Trilogy200, universal porting block	
Reusable circuits		
1052064	Trilogy active w/PAP exhalation circuit, adult	A4618
1052065	Trilogy active w/PAP exhalation circuit, ped	A4618
1052066	Trilogy passive exhalation circuit, adult	A4618
1052127	Trilogy passive exhalation circuit, pediatric	A4618
1068837	Trilogy200 active w/Flow sensor circuit, adult	A4618
1068838	Trilogy200 active w/Flow sensor circuit, pediatric	A4618



### MPV circuit support system

An interlocking support system combines with Trilogy's breathing circuit to make access easy for patients. The support system can be attached to various surfaces in many configurations with the use of a clamp designed for use on flat or rounded surfaces. The interlocking support circuit can be adjusted to enhance comfort and accessibility.



Under Medicare, all parts and accessories are included in the monthly rental of the device and are not separately billable. Supplies may be separately billable for patient-owned ventilators. Consult the payer policy manual for specific details.



Part number		HCPCS
<b>Adult circuits</b>		
1069467	Adult passive circuit with fixed exhalation port, 10 pk	A4618
1069468	Adult passive circuit with fixed exhalation port, 20 pk	A4618
1073228	Adult passive circuit with swivel design exhalation port, 10 pk	A4618
1073236	Adult passive circuit with swivel design exhalation port, 20 pk	A4618
1069443	Adult passive circuit with water trap, 10 pk	A4618
1069444	Adult passive circuit with water trap, 20 pk	A4618
1073227	Adult passive circuit, heated wire, 10 pk	A4618
1073235	Adult passive circuit, heated wire, 20 pk	A4618
1073223	Adult passive circuit, no water trap, 10 pk	A4618
1073232	Adult passive circuit, no water trap, 20 pk	A4618
1073222	Adult active circuit, heated wire, 10 pk	A4618
1073231	Adult active circuit, heated wire, 20 pk	A4618
1102866	MPV circuit, 10 pk	A4618
<b>Pediatric circuits</b>		
1073224	Pediatric passive circuit with water trap, 10 pk	A4618
1073233	Pediatric passive circuit with water trap, 20 pk	A4618
1073225	Pediatric passive circuit, heated wire, 10 pk	A4618
1073234	Pediatric passive circuit, heated wire, 20 pk	A4618
1073220	Pediatric active circuit with water trap, 10 pk	A4618
1073229	Pediatric active circuit with water trap, 20 pk	A4618
1073221	Pediatric active circuit, heated wire, 10 pk	A4618
1073230	Pediatric active circuit, heated wire, 20 pk	A4618
<b>Circuit adapters</b>		
1073902	Flexible trach adapter with 22 mm cuff (connects to Whisper Swivel II and swivel and fixed disposable exhalation ports), 10 pk	
1073904	Flexible trach adapter with 22 mm cuff (connects to Whisper Swivel II and swivel and fixed disposable exhalation ports), 20 pk	
1073863	Flexible trach adapter with 15 mm cuff (connects to the flow sensor or the active exhalation valve), 10 pk	
1073903	Flexible trach adapter with 15 mm cuff (connects to the flow sensor or the active exhalation valve), 20 pk	
1065659	Active exhalation valve with PAP assembly	
1061398	Trilogy 200, universal active PAP tube adapter	
1050408	Trilogy 200, proximal flow sensor	
1050534	Trilogy 200 Proximal flow sensor, 10 pk	
332113	Whisper Swivel II	
312149	Disposable fixed exhalation port, 10 pk	
1084428	Disposable swivel passive exhalation port, 10 pk	
1084429	Disposable swivel passive exhalation port, 20 pk	

## Trilogy100 and Trilogy200

Part number		HCPCS
<b>Carry cases</b>		
1033709	Trilogy travel carrying case	
1040420	Trilogy carrying case, in-use	
<b>Battery packs</b>		
1072182	Detachable battery charger	A4613
1043570	Trilogy detachable battery pack, dom/intl	A4611
<b>Filters</b>		
1029330	Trilogy filter, 2 pk	
C06417	Filter, bacteria, 22 mm, single use	
C06418	Filter, bacteria, 22 mm, single use, 50 pk	
<b>Stands</b>		
1047410	Trilogy small products roll stand	
1094043	Trilogy circuit support arm	
1074004	RP-mount plate, Trilogy roll stand	
1074137	RP screw, MNT plate, Trilogy roll stand, 4 pk	
<b>Humidifiers</b>		
22135	Humidifier, Fisher & Paykel, MR410	E0550
1074173	Fisher & Paykel HC500 humidifier, domestic	E0550
1074174	Fisher & Paykel airway temperature probe	
1074175	Fisher & Paykel heated hose humidifier adapter	
<b>Pulse oximetry</b>		
1075383	Masimo high performance/low power pulse oximeter	
1080964	Trilogy interface kit: interface adapter, SD card and mailer, carry case	
1081984	Masimo LNCS DCI, adult reusable sensor	
1081515	Masimo LNCS DCIP, pediatric reusable sensor	
1081517	Masimo LNCS ADTX, adult disposable sensor, 20 bx	
1081518	Masimo LNCS PDTX, pediatric disposable sensor, 20 bx	
1081985	Masimo LNCSYI, multisite reusable sensor	
1081516	Masimo LNCS DBI, adult soft reusable sensor	
1082474	Masimo LNC-4-Ext, 4 ft extension cable	
1019531	Standard multisite wrap, 100 pk	
1019532	Clean shield multisite wrap, Masimo, 100 pk	
<b>Accessories</b>		
1102862	MPV support system	
1040390	Trilogy, O <sub>2</sub> quick-disconnect insert, 2 pk	
1021671	Test lung	
1045587	Inlet airports assembly	

Trilogy accessories are separately billable only when the ventilator is patient owned.



# Airway clearance

For use by caregivers and patients, in the hospital or in the home, CoughAssist T70 offers a gentle and more comfortable alternative to traditional suctioning methods.

Not all products are available in all countries. Check with Philips Respironics customer service or your local distributor.





# Clearly a better solution

**CoughAssist** is an effective, more comfortable way to clear airways

CoughAssist T70 is a noninvasive therapy that safely and consistently removes secretions in patients with an ineffective ability to cough. It clears secretions by gradually applying a positive pressure (insufflation) to the airway and then rapidly shifting to negative air pressure (exsufflation). This in-exsufflation shift in pressure creates a high expiratory flow that simulates a deep, natural cough. Instead of introducing a suction catheter into the airway, air is delivered noninvasively through a facemask, a mouthpiece, or through an adapter that allows the device to function with an endotracheal or tracheostomy tube.

#### Key features

- Digital display and push button controls
- AC, DC and optional battery power options
- Cough-Trak automatic sensitivity
- DirectView data management

To learn more about CoughAssist T70 visit [www.philips.us/coughassistt70](http://www.philips.us/coughassistt70)



Part number		HCPSC
1098160	CoughAssist T70	E0482
1098159	CA E70 (Canada)	E0482
Circuits		
1090830	CA 70 series circuit, infant, 6 ft	A7020
1090831	CA 70 series circuit, toddler, 6 ft	A7020
1090832	CA 70 series circuit, adult small, 6 ft	A7020
1090833	CA 70 series circuit, adult medium, 6 ft	A7020
1090834	CA 70 series circuit, adult large, 6 ft	A7020
1090835	CA 70 series circuit, trach, 6 ft	A7020
1090836	CA 70 series circuit, mouthpiece, 6 ft	A7020
1098403	CA 70 series circuit, infant, 9 ft	A7020
1098404	CA 70 series circuit, toddler, 9 ft	A7020
1098405	CA 70 series circuit, adult small, 9 ft	A7020
1098407	CA 70 series circuit, adult medium, 9 ft	A7020
1098408	CA 70 series circuit, adult large, 9 ft	A7020
1098409	CA 70 series circuit, trach, 9 ft	A7020
1098410	CA 70 series circuit, mouthpiece, 9 ft	A7020
Masks		
1046861	CoughAssist, facemask, adult large	
1044187	CoughAssist, facemask, adult medium	
1044188	CoughAssist, facemask, adult small	
1044190	CoughAssist, facemask, infant	
1044189	CoughAssist, facemask, toddler	
Miscellaneous		
1046860	CoughAssist, disposable filter	
1099533	CoughAssist, mouthpiece	
1098720	CA 70 series, water trap	
1098655	CA 70 series, roll stand	
1098656	CA 70 series, roll stand bracket	
1099034	CA T70 oximetry interface kit	
1098718	CA E70 oximetry interface kit	
1059017	CA 70 series, foot pedal	
1098884	CA 70 series, carrying case	
1051801	SD card and case	
1099035	CA 70 series, circuit retainer	
1032584	Pollen filter, single pack	
1043570	Trilogy detachable battery pack, dom/intl	
1099036	CA 70 series, battery cover	
1072182	Detachable battery charger	
1099427	RP-power cord, 8 ft, C17, NEMA 1-15P	
1067424	Trilogy/CoughAssist automotive DC power adapter	



Accessories are not separately billable to the Medicare program, with the exception of replacement circuits. Consult your local Medicare DME MAC local coverage determination (LCD) for details. Commercial insurance plans may provide separate coverage for accessories. Consult insurance company policy manual for more information.



# Respiratory drug delivery

Respiratory Drug Delivery provides innovative drug delivery technologies for the treatment of respiratory and non-respiratory diseases. Our easy-to-use products include compressor nebulizer systems, holding chambers, and peak flow meters, as well as intelligent drug delivery systems that enable the delivery of precise and reproducible doses of drug.

Not all products are available in all countries. Check with Philips Respironics customer service or your local distributor.

# Compressor nebulizers

## InnoSpire Essence compressor nebulizer system

Compact, lightweight and economical, the InnoSpire Essence compressor nebulizer system combines with proven SideStream technology to provide fast and efficient aerosol drug delivery.

To learn more visit  
[www.philips.us/essence](http://www.philips.us/essence)

Part number		HCPCS
1099966	InnoSpire Essence with SideStream disposable nebulizer	E0570
1100312	InnoSpire Essence with SideStream disposable and reusable nebulizers	E0570
1100313	InnoSpire Essence with SideStream disposable and reusable nebulizers, Tucker the Turtle mask, and carrying case	E0570



## InnoSpire Elegance compressor nebulizer system

The InnoSpire Elegance compressor combined with the SideStream nebulizer provides a powerful, high-quality compressor nebulizer system that is durable, easy to use, and easy to maintain. It provides the freedom to move, helping your patients build solid compliance behaviors – in a size that easily fits into their lives.

To learn more visit  
[www.philips.us/elegance](http://www.philips.us/elegance)

Part number		HCPCS
1100134	InnoSpire Elegance with SideStream disposable nebulizer	E0570
1099969	InnoSpire Elegance with SideStream disposable and reusable nebulizers	E0570
1102828	InnoSpire Elegance with SideStream Plus	E0570

Accessories		HCPCS
<b>Filters</b>		
1102088	Replacement filter, 4 per pk	A7013
<b>Miscellaneous</b>		
1026568	Premium carrying case	
1026641	Standard carrying case	
1101452	InnoSpire carrying case	



# Compressor nebulizers

## InnoSpire Deluxe compressor nebulizer system

For use with aerosol medications, InnoSpire Deluxe incorporates built-in conveniences that encourage adherence. It includes a convenient storage compartment to hold extra medicine or nebulizer supplies, the attached power cord tucks neatly into the storage area, and a sturdy handle makes it easy to carry. InnoSpire Deluxe is designed with compliance in mind.

Part number		HCPCS
1110063	InnoSpire Deluxe compressor with SideStream disposable nebulizer (6/cs)	E0570
1102088	Filter, 4 pk	A7013



## Sami the Seal compressor nebulizer system

The Sami the Seal pediatric compressor is shaped like a seal to appeal to children, their parents, and pediatricians, and to support aerosol therapy compliance. The compressor partners with highly efficient SideStream nebulizers to provide fast and friendly treatment to pediatric patients in the home.

Part number		HCPCS
1082456	Sami compressor with SideStream disposable nebulizer and Tucker the Turtle mask	E0570
1082731	Sami compressor with SideStream disposable and reusable nebulizers and Tucker the Turtle mask	E0570

Accessories		HCPCS
1083128	Replacement filters, 5 pk	A7013
1099821	Carrying case	





# Compressor nebulizers

## MicroElite compressor nebulizer system

Proven jet nebulizer technology and a small, compact size combine to make MicroElite uniquely portable and convenient. Half the size of other common portables, it is easy to use and transport. It includes a DC adapter and an optional rechargeable lithium polymer battery. Coupled with quick treatment times, MicroElite provides nebulizer patients with portability and speed that combine to optimize treatment and promote active lifestyles.

Part number		HCPCS
RDD490	MicroElite compressor nebulizer system, standard	E0570
RDD491	MicroElite compressor nebulizer system, deluxe	E0570

Accessories		HCPCS
<b>Miscellaneous</b>		
RDD499	MicroElite rechargeable lithium polymer battery	
1041431	MicroPlus high-efficiency reusable nebulizer	A7005
1041539	MicroElite AC power adapter	
1041540	MicroElite 12VDC (car) adapter	
1041534	MicroElite replacement air filter packs, 4 per pk	A7013
1049792	MicroElite pediatric mask	A7015
1049793	MicroElite adult mask	A7015



# Compressor nebulizers

## MiniElite compressor nebulizer system

The MiniElite compressor nebulizer system is compact, lightweight, and designed for patient convenience. It can be powered by three versatile power sources – 110VAC, 12VDC (car), and an optional rechargeable lithium ion battery. The freedom to choose when and where patients take their aerosol treatments makes the MiniElite an ideal daily traveling companion.

Part number		HCPCS
RDD482	MiniElite compressor nebulizer, standard	E0570
RDD483	MiniElite compressor nebulizer with battery	E0570

Accessories		HCPCS
Miscellaneous		
RDD489	MiniElite rechargeable lithium ion battery	
1036899	MiniElite AC power adapter	
1036519	MiniElite 12VDC (car) adapter	
1036898	MiniElite replacement air filter packs, 4 pk	A7013



# Nebulizers

## SideStream disposable nebulizer

The SideStream disposable nebulizer offers fast, efficient, and economical drug delivery. It provides consistent drug delivery time after time and can be used with all commonly prescribed respiratory medications. It can be used for two weeks before replacing.

Part number		HCPCS
HS800	SideStream disposable nebulizer, 50 pk	A7003
1028219	SideStream disposable cup, 50 pk	A7004
1036160	SideStream disposable with adult mask, 50 pk	A7003,A7015
1041906	SideStream disposable with pediatric mask, 50 pk	A7003,A7015



## SideStream reusable nebulizer

The SideStream reusable nebulizer provides maximum aerosol quality and respirable output for the delivery of most commonly prescribed medications. The unique venturi design offers faster drug delivery and shorter treatment times. It can be used for six months before replacing.

Fast treatment times – delivers 2.5 ml of drug output within seven minutes.

Highly efficient – 80 percent of drug output volume is below five microns in size.

Part number		HCPCS
HS860	SideStream reusable nebulizer, 10 pk	A7005
HS860-050	SideStream reusable nebulizer, 50 pk	A7005



## SideStream Plus reusable nebulizer

The SideStream Plus breath-enhanced nebulizer features an easy-action inspiratory valve that opens on inspiration to boost medication delivery and closes on exhalation, preserving medication and reducing waste. Significantly less pressure is required to open the inspiratory valve compared to other breath-enhanced nebulizers, making treatments easier for pediatric and elderly patients. And its proven SideStream five-jet venturi design provides a high respirable fraction and shortens treatment times.

Part number		HCPCS
HS870	SideStream Plus breath-enhanced nebulizer, 10 pk	A7005



# Nebulizer kits

## SideStream custom nebulizer kit

SideStream provides superior performance and fast treatment times. SideStream utilizes five aerosol-producing jet holes compared to a conventional nebulizer's one. The SideStream kit provides a two-month supply of product in every package.

Part number		HCPCS
1028125	SideStream custom nebulizer kits, 25 pk	A7003 and A7004

*Includes two complete SideStream disposable nebulizers and two replacement nebulizer cups.*





# Nebulizer masks

## SideStream masks

The SideStream mask is recommended for use with the SideStream high-efficiency nebulizers. The reusable mask is contoured to fit the face. The pediatric mask is recommended for ages five and under; the adult size for ages six and up.

Part number		HCPCS
1025529	SideStream pediatric mask, 50 pk	A7015
1025531	SideStream adult mask, 50 pk	A7015
1052521	SideStream Plus adult mask	A7015



## SideStream pediatric eye-safe masks

Sami the Seal and Tucker the Turtle child-friendly masks are designed to appeal to children and have been shown to reduce aerosol deposition in and around the eyes. Used with SideStream nebulizers, Sami and Tucker pediatric masks are made of soft, flexible material for comfort and fit.

Part number		HCPCS
HS880	Sami the Seal pediatric mask	A7015
1044382	Tucker the Turtle pediatric mask	A7015





## A cut above the rest

### OptiChamber Diamond antistatic valved holding chamber

#### **The "one for all" chamber**

In the hospital or at home, its uniquely designed stepped mouthpiece and detachable masks make the OptiChamber Diamond suitable for patients of all ages.

#### **The chamber that fits into your patients' lives**

OptiChamber Diamond is smaller than most conventional VHCs and is designed for patient comfort. Its stepped mouthpiece facilitates transfer from pediatric mask to mouthpiece and its flat bottom maintains stability when the VHC is not in use.

#### **Key features**

- Intuitive, ergonomic design
- Low resistance expiratory valve
- Flexible pMDI adapter
- Detachable LiteTouch facemasks with soft-seal technology

Part number		HCPSC
1077478	OptiChamber Diamond, 10 pk	A4627
1079823	OptiChamber Diamond with LiteTouch small mask, 10 pk	A4627
1079826	OptiChamber Diamond with LiteTouch medium mask, 10 pk	A4627
1079827	OptiChamber Diamond with LiteTouch large mask, 10 pk	A4627
1098213	OptiChamber Diamond, 10 pk, Canada	
1098214	OptiChamber Diamond with Lite Touch small mask, 10 pk, Canada	
1098215	OptiChamber Diamond with LiteTouch medium mask, 10 pk, Canada	
1098216	OptiChamber Diamond with LiteTouch large mask, 10 pk, Canada	

### OptiChamber Diamond antistatic valved holding chamber (VHC)


The OptiChamber Diamond antistatic valved holding chamber (VHC) is designed for ease of use, to help maximize delivery of pMDI medications to the lungs, and encourage compliance for patients of all ages, at home or in the hospital. Diamond is one of the smallest full-function VHCs available. Its low-resistance valves open freely, even under low pediatric pressures and flows, its flat bottom maintains stability when not in use, its clear chamber allows visualization of medication, and its expiratory valve allows caregivers to count patient breaths. Diamond's stepped mouthpiece better accommodates smaller mouths and helps children easily graduate from mask to mouthpiece. Paired with the comfort and distinctive soft-seal interface of the LiteTouch detachable mask (available in multiple sizes), Diamond promotes aerosol therapy comfort and compliance.

To learn more visit [www.philips.com/optichamberdiamond](http://www.philips.com/optichamberdiamond)


**OptiChamber Diamond is designed for your patients, including those who have:**

- **Poor coordination**  
The antistatic material suspends the aerosol longer, meaning your patient has more time to inhale their medication<sup>[1]</sup>.
- **A low inspiratory flow rate**  
Low resistance valves let the patient breath easily through the chamber with minimal resistance.
- **Use a mask**  
Simply attach the LiteTouch facemask, which has soft-seal technology and requires minimal pressure to maintain a proper seal<sup>[2]</sup>.
- **Their medication administered by someone else**  
The exhalation valve visually confirms exhalation, helping the caregiver to count patient breaths and coordinate actuation.

## New, smaller and antistatic



### LiteTouch VHC face mask



1. Mitchell, J., et al. "Electrostatics and Inhaled Medications: Influence on delivery via pressurized metered-dose inhalers and add-on devices," *Resp. Care*. March 2007, Vol 52 No 3.

2. Hsu, W., et al. (2011), "Evaluation of delivery efficiency from valved holding chambers with facemasks," *Resp. Drug. Deliv. Europe* 2011, 2, pp 505-510. Respironics, OptiChamber and LiteTouch are trademarks of Koninklijke Philips Electronics N.V. and its affiliates. All rights reserved.

# Spacers, chambers, and asthma kits

## ProChamber valved holding chamber

ProChamber is the professional model valved holding chamber for the acute respiratory care setting. Its proven clinical performance and great value make ProChamber the choice for the budget-conscious clinician.

Part number		HCPCS
HS2003-050	ProChamber valved holding chamber, 50 pk	A4627



## AsthmaPACK

The AsthmaPACK adult and children's personal asthma care kit for patient adherence implements three key elements of the NAEPP asthma treatment guidelines with one prescription-peak flow monitoring, MDI drug delivery, and asthma education. One prescription will give your patients a comprehensive, clinician-directed, patient-managed asthma plan with the tools needed to better understand their asthma, monitor their conditions objectively, and take their prescribed medications correctly and consistently.

Part number	
1099014	Adult asthma care kit, 10 pk
1099015	Children's asthma care kit, 10 pk

*Includes a peak flow meter, spacer, instructional video, and educational self-help booklet.*



## AsthmaPACK children's asthma care kit

AsthmaPACK children's asthma care kit makes learning to manage asthma child's play. This just-for-kids kit makes learning fun with engaging cartoon characters and age-appropriate devices. An exciting feature of this kit is the interactive bilingual (English/Spanish) coloring and activity storybook, "Tucker Tackles Asthma."





# Chambers and masks

## LiteTouch valved holding chamber mask

LiteTouch fuses a hard shell to an exclusive soft seal that attaches onto the mouthpiece of the OptiChamber Diamond valved holding chamber (VHC). Gently touching the mask to the face creates a cushion-like seal that quickly contours to the face, making it easy to fit, especially for small children. Coupled with the Diamond VHC, LiteTouch's soft-seal technology minimizes leaks and helps to deliver optimum aerosol therapy for asthma and other respiratory illnesses.

Part number	HCPCS
1082714	LiteTouch mask small, 10 pk
1082715	LiteTouch mask medium, 10 pk
1082716	LiteTouch mask, large, 10 pk



## Silicon masks

Our detachable valved holding chamber masks provide an effective seal and conform to fit different facial geometries. Built-in reinforcements help prevent the mask from collapsing and assist in maintaining a seal, especially when your pediatric patients are not cooperative. Available in small (0 – 18 months), medium (1 – 6 years), and large (6 years +) sizes.

Part number	HCPCS
HS81111-010	OptiChamber small face mask, 10 pk
HS81211-010	OptiChamber medium face mask, 10 pk
HS81311-010	OptiChamber large face mask, 10 pk



# Peak flow meters

## ASSESS

Measure airway obstruction easily, accurately, and cost-effectively. Peak expiratory flow rate (PEFR) provides an objective, clinically-relevant measurement of airflow – at a fraction of the cost, bulk, and inconvenience of conventional spirometry. And, when seconds count, you can count on ASSESS to deliver those measurements with superior accuracy and reproducibility. Everything about ASSESS, from its easy-to-use design, to its easy-to-read scale, to its transparent construction that encourages regular cleaning, to its virtually indestructible polycarbonate body, is designed to make it work better and last longer.

Part number	HCPCS
HS710-012	A4614

ASSESS full range peak flow meter, 12 pk

Unit includes an adult and pediatric mouthpiece.



## AsthmaCheck

AsthmaCheck's easy-grip handle directs hands away from the airflow path and ensures accurate measurements every time. Like all of our peak flow meters, it has a built-in three-zone management system and an easy-to-read scale which is suitable to all patients (60 – 810 l/min).

Part number	HCPCS
HS740-012	A4614

AsthmaCheck, 12 pk



# Peak flow meters

## AsthmaMentor

Turn our peak flow meter into your asthma action plan and teach your patients better asthma management in minutes. AsthmaMentor has been designed to help clinicians teach patients what to do when their asthma conditions change. AsthmaMentor's unique features make the most of the time you have with your patients. The exclusive built-in AutoZone system automatically calculates green, yellow, and red zones when the top of the color bar is set at a patient's personal best or predicted average PEF. The asthma action plan sticker, which provides color-coded areas where you can write medication instructions, can be affixed to the back of the meter. Your instructions are then easily matched with daily peak flow readings, helping your patient comply with them.

Part number		HCPCS
HS742-010	AsthmaMentor, 10 pk	A4614



## Personal Best peak flow meter

Wherever Personal Best goes – in the hospital, office, or home – it sets a new standard for convenience, compliance, and confidence. And Personal Best goes anywhere...easily. A lightweight, self-contained, portable design puts peak flow monitoring at your fingertips – or your patients' – all day long. Plus, a rugged, sterilizable construction helps it stand up to the roughest use day after day. Personal Best encourages compliance with its comfortable oval mouthpiece and easy-to-use built-in handle. And it improves patient communication with its high-legibility scale and NAEPP-endorsed three-zone management system. Available in both full range (60 – 810 l/min) and low range (50 – 390 l/min).

Part number		HCPCS
HS755-012	Personal Best full range peak flow meter, 12 pk	A4614
HS756-012	Personal Best low range peak flow meter, 12 pk	A4614



# Mouthpieces

## Disposable and one-way valve safety mouthpieces

Safety mouthpieces reduce the risk of cross-contamination during multi-patient screening with Personal Best, ASSESS, AsthmaMentor, and AsthmaCheck peak flow meters. The one-way valve prevents patients from inadvertently breathing in air through their meters during PEFr testing.

Part number	
HS711-100	Disposable adult mouthpieces for ASSESS, 100 pk
HS712-100	Disposable pediatric mouthpieces for ASSESS, 100 pk
HS715-100	Disposable pediatric plastic mouthpieces for all Respironics' peak flow meters, 100 pk
HS757-100	Disposable mouthpieces for Personal Best, AsthmaMentor, and AsthmaCheck, 100 pk
HS713-200	One-way valve safety mouthpieces for ASSESS, 200 pk
HS714-200	One-way valve safety mouthpieces for Personal Best, AsthmaMentor, and AsthmaCheck, 200 pk




# Accessories

## OptiVent MDI in-line spacer

The OptiVent MDI in-line spacer is designed to optimize drug delivery and reduce fluid accumulation. Convenient and easy to use, it does not need to be expanded when used with a ventilator circuit.

Part number	
HS762-010	OptiVent, 10 pk
HS763-020	OptiVent with 15mm adapters, 20 pk





# Inspiratory muscle trainers

## Pflex IMT

Pflex is a highly-effective inspiratory muscle trainer designed to increase respiratory muscle strength. The Pflex trainer is based on the principle of resistance to intake of air into the lungs which makes respiratory muscles work harder. Resistance is controlled by six fixed orifice settings.

Part number		HCPCS
HS553-010	Pflex IMT, 10 pk	E1399



## Threshold IMT

Provides consistent and specific pressure for inspiratory muscle strength and endurance training, regardless of how quickly or slowly patients breathe. This device incorporates a flow-independent, one-way valve to ensure consistent resistance and features an adjustable specific pressure setting (in cm H<sub>2</sub>O) to be set by a healthcare professional. When patients inhale through Threshold IMT, a spring-loaded valve provides resistance that exercises respiratory muscles through conditioning.

Part number		HCPCS
HS730-010	Threshold IMT, 10 pk	E1399



## Threshold PEP

Threshold PEP is used for airway clearance, bronchial hygiene, or as an alternative to chest physical therapy. Threshold PEP incorporates a flow-independent, one-way valve to ensure consistent resistance with adjustable specific pressure settings (in cm H<sub>2</sub>O) which are set by the healthcare professional. When patients exhale through Threshold PEP, the resistive load creates positive pressure that helps open the airways and allows mucus to be expelled during "huff" coughing (forced expiratory technique).

Part number		HCPCS
HS735-010	Threshold PEP, 10 pk	E1399





# Oxygen

Philips Respironics Right Fit oxygen therapy products are intended to give patients the freedom they want while helping homecare providers reduce expenses and make the most of shrinking reimbursements.

Not all products are available in all countries. Check with Philips Respironics customer service or your local distributor.

# Concentrators

Part number		HCPCS
1020000	EverFlo oxygen concentrator	E1390
1020001	EverFlo oxygen concentrator with oxygen percentage indicator (OPI)	E1390
1039362	EverFlo oxygen concentrator, transfill	E1390
1039363	EverFlo oxygen concentrator with OPI, transfill	E1390

Units include a patient manual.



## EverFlo concentrator

Delivers a great performance and demands so little in return. It's small, quiet, and designed to be easier and less expensive to operate. Plus, EverFlo reduces scheduled maintenance visits. So now you can control oxygen therapy costs and deliver all the features you and your patients want with one 31-pound concentrator that's bound to be the star of your business.

Accessories	
<b>Filters</b>	
1038831	Compressor intake filter
520	Micro disk filter
<b>Tubing</b>	
1039642	Humidifier connector tube
<b>Miscellaneous</b>	
H537-1	Stand-alone pediatric flow meter
H644	Integrated pediatric flow meter
528	Flow meter locking kit
H649	Caster brake kit
1038571	Warranty, two-year with compressor

Part number		HCPCS
1020015	EverFlo Q oxygen concentrator	E1390
1020014	EverFlo Q oxygen concentrator with OPI	E1390
1039365	EverFlo Q oxygen concentrator, transfill	E1390
1039364	EverFlo Q oxygen concentrator with OPI, transfill	E1390

Units include a patient manual.



## EverFlo Q ultra quiet 40 decibel stationary oxygen concentrator

Our compact, 31-pound EverFlo Q is the first ultra-quiet concentrator available at an affordable price. The new EverFlo Q combines the small size, light weight, and low maintenance of our popular EverFlo oxygen concentrator, with the added value of an ultra-quiet sound level. With EverFlo Q, you get the benefits of an ultra-quiet concentrator without the burden of paying an ultra-large premium. And, that's good for you and your patients.

Accessories	
<b>Filters</b>	
1038831	Compressor intake filter
520	Micro disk filter
<b>Tubing</b>	
1039642	Humidifier connector tube
<b>Miscellaneous</b>	
H537-1	Stand-alone pediatric flow meter
H644	Integrated pediatric flow meter
528	Flow meter locking kit
H649	Caster brake kit
1038571	Warranty, two-year with compressor



# Concentrators



## Life made simple

### Philips Respironics SimplyFlo stationary oxygen concentrator

For nocturnal patients who may be refusing or resisting therapy because of the size, weight and lifestyle disruptions that can accompany the use of traditional concentrators, SimplyFlo weighs less than nine pounds and is designed to enhance compliance. SimplyFlo is the smallest and lightest continuous flow stationary oxygen concentrator available. It operates on external power only and has no internal battery. An optional travel accessory case allows oxygen patients to easily transport SimplyFlo between locations.

#### Key features

- Weighs 8.5 pounds
- Provides .5 to 2 LPM continuous flow output
- Operates on AC power (no internal battery)
- Uses less than 1/3 the electricity of traditional stationary concentrators
- Reduces many lifestyle disruptions related to larger stationary concentrators
- Simplifies transport for travel or within the home
- Helps nocturnal oxygen users maintain an active lifestyle during the day

Part number		HCPCS
1076773	SimplyFlo oxygen concentrator (includes concentrator, AC power cord and user manual	E1390

Accessories		
1076830	SimplyFlo travel case	
1082664	AC power cord	
1082661	AC power supply	
1076809	User manual	
1076829	Humidifier bottle holder	
1101172	Humidifier connecting tube kit	
Warranties		
1097623	SimplyFlo two-year total warranty	
1097577	SimplyFlo three-year total warranty	



## Another solution for nocturnal patients

The optional travel accessory case allows patients to easily transport SimplyFlo between locations.



Once again, Philips Respironics is bringing innovative solutions to your business and your patients. With SimplyFlo, we've created the world's smallest and lightest stationary oxygen concentrator.

SimplyFlo represents a completely new option for nocturnal patients who may be refusing or resisting therapy because of the size, weight and lifestyle disruptions that can accompany the use of traditional concentrators.

### Designed to enhance compliance

For many patients new to therapy who only require nocturnal treatment, the intimidating size, weight, and sound of a traditional oxygen concentrator can be a large factor in their rejection of therapy.

Until SimplyFlo, there has not been a significantly lighter, more transportable stationary concentrator designed to help nocturnal oxygen patients improve the adoption of therapy and enhance compliance.



# Concentrators

Part number		HCPCS
M10600	Millennium 10-liter oxygen M10 concentrator	E1390
M10600-80	Millennium 10-liter oxygen M10 concentrator, Canada	
M10605	Millennium 10-liter oxygen M10 concentrator with OPI, Canada	E1390
M10605-80	Millennium 10-liter oxygen M10 concentrator with OPI, Canada	
M10600-10	Millennium 10-liter oxygen concentrator, transfill	E1390
M10605-10	Millennium 10-liter oxygen concentrator with OPI, transfill	E1390

Both units include a patient manual.

Accessories		HCPCS
1014653	High flow starter kit <i>Includes cannula, 25' tubing, connector</i>	A4615*, A4616*
1014654	High flow nasal cannula	A4615*
1017796	Inlet filter replacement kit	
1017768	High flow oxygen tubing, 25'	A4616*
1014655	Oxygen tubing connector	

\*Accessories are not separately billable to the Medicare program, with the exception of replacement circuits. Consult your local Medicare DME MAC local coverage determination (LCD) for details. Commercial insurance plans may provide separate coverage for accessories. Consult insurance company policy manual for more information.



## Millennium M10 concentrator

Delivers up to 10 LPM of oxygen reducing the delivery costs associated with 5 LPM and above oxygen patients. The M10 is engineered to reduce the cost of providing oxygen at higher liter flow and has fewer parts than any other concentrator on the market. Fewer parts, combined with greater reliability, help you to reduce expenditures on accessories and unnecessary inventory.

# Portable concentrators



## Simplify and SimplyGo

Now one POC meets the portable oxygen needs of nearly all of your products.

SimplyGo is Philips first oxygen concentrator that offers the benefits of a continuous flow mode in a lightweight, easy-to-carry design of 4.5 kg. No need to compromise between an effective oxygen therapy output for your patients and the usual weight of the device. SimplyGo offers a solution that fulfills the needs of most of your oxygen-therapy patients. With SimplyGo, patients have the flexibility to be away for the day, for work, for a weekend or a long holiday.

To learn more visit [www.simplygo.respironics.com](http://www.simplygo.respironics.com)

### Key features

- Attractive, sleek design
- Easy-to-remove battery
- Intuitive, easy-to-read screen
- Convenient carrying case
- Handy fold-up cart
- Detachable accessory bag

Part number		HCPCS
1068987	SimplyGo portable oxygen concentrator	E1390, E1392
1068787, 1020000	SimplyGo portable oxygen concentrator with EverFlo non-OPI oxygen concentrator	E1390, E1392
1068787, 1020001	EverFlo OPI oxygen concentrator	E1390, E1392
1068787, 1020015	EverFlo Q non-OPI oxygen concentrator	E1390, E1392
1068787, 1020014	SimplyGo portable oxygen concentrator with EverFlo Q with OPI oxygen concentrator	E1390, E1392

Accessories		
1082663	SimplyGo carry case	
1083692	SimplyGo DC power supply	
1082661	SimplyGo AC power supply and cord	
1083693	SimplyGo airline power cord	
1083696	SimplyGo accessory bag	
1106598	SimplyGo shoulder strap	
1083699	SimplyGo humidifier pouch	
1101172	SimplyGo humidifier connecting tube kit	
1101602	SimplyGo humidifier connecting tube and pouch	
1082662	SimplyGo lithium ion battery with handle	
900-103	Battery charger	
1074885	Mobile cart	
1083700	User manual	
Warranties		
1097623	SimplyGo two-year total warranty	
1097577	SimplyGo three-year total warranty	
1097778	SimplyGo five-year total warranty	



### SimplyGo portable oxygen concentrator

SimplyGo is the only portable oxygen concentrator (POC) to offer continuous flow and pulse-dose delivery in a single device weighing only 10 pounds. That means you can manage the therapy and lifestyle needs of nearly all of your patients with just one POC. Designed to be rugged and easy to use, SimplyGo is truly a new breed of lightweight and attractive portable concentrator.

## SimplyGo at a glance

**Simple-to-select oxygen delivery mode**

**Simple-to-use touch-pad interface**

**Lightweight - 4.5 kg**

The only lightweight device able to deliver both, pulse and continuous oxygen flow.

**Simple-to-carry anywhere, anytime!**

**Simple-to-replace/remove rechargeable battery.**

See product specifications table for more information about battery life.



The United States Federal Aviation Administration (FAA) has approved the SimplyGo for in-flight use on commercial airliners.



# Home oxygen filling system



Part number		HCPCS
1057100	UltraFill station 120V	K0738
1057101	UltraFill system with 2 cylinders, 2,000 psi and a cylinder bag, 120V, Canada	K0738
1074932	UltraFill system with 2 cylinders, 2,000 psi and a cylinder bag, 120V <i>3,000 psi cylinders available for additional charge.</i>	K0738
1039362	EverFlo 120V, transfill	E1390
1039363	EverFlo 120V OPI, transfill	E1390
1039365	EverFlo Q 120V, transfill	E1390
1039364	EverFlo Q 120V OPI, transfill	E1390
M10600-10	Millennium 120V, transfill	E1390 (with QF or QG modifier)
M10605-10	Millennium 120V OPI, transfill	E1390 (with QF or QG modifier)

## Accessories

### Cylinders

1065705	M4 cylinder, US, 2000 psi
1065706	M6 cylinder, US, 2000 psi
1065707	M9 cylinder, US, 2000 psi
1065708	MD15 cylinder, US, 2000 psi
1065709	ME24 cylinder, US, 2000 psi
1065710	M4 cylinder, Canada, 2000 psi
1065711	M6 cylinder, Canada, 2000 psi
1065712	M9 cylinder, Canada, 2000 psi
1065713	MD15 cylinder, Canada, 2000 psi
1065714	ME24 cylinder, Canada, 2000 psi
1065715	MB08 cylinder, US, 3000 psi
1065716	MC13 cylinder, US, 3000 psi
1065717	MD22 cylinder, US, 3000 psi
1065718	ME36 cylinder, US, 3000 psi
1065719	MB08 cylinder, Canada, 3000 psi
1065720	MC13 cylinder, Canada, 3000 psi
1065721	MD22 cylinder, Canada, 3000 psi
1065722	ME36 cylinder, Canada, 3000 psi

### Accessories

1065697	Chad Therapeutics, Bonsai oxygen conserving device
1065698	Chad Therapeutics, 3100 pediatric regulator
1065699	Chad Therapeutics, 4800 regulator
1065700	Chad Therapeutics, 3100 regulator
1065701	Small cylinder bag (M9/MC13)
1065702	Medium cylinder bag (M4/M6/MB08)
1065703	Large cylinder bag (MD15/MD22)
1065724	Cylinder cart (MD15/ME24/MD22/ME36)
1065696	UltraFill transfer tubing

## UltraFill home oxygen filling system

Is an advanced home oxygen system that combines a stationary oxygen concentrator, filling station, and high-capacity cylinders to meet the needs of a wide range of oxygen patients, including those who are highly active or require continuous flow oxygen. UltraFill also has the needs of your business in mind with features intended to save you money by using existing inventory and lowering your operating costs.

# Oximeters

## Part number

920	Hand-held oximeter
920M	Hand-held oximeter with memory

Both include a 936 adult finger clip sensor, six AA batteries, and a manual.

## Accessories

### Hand-held oximeter

100-0433-00	Battery door
927	Extension cable, 3'
1008130	Extension cable, 10'
924	Manual
1007332	Carrying case

### Hand-held oximeter with memory

H920M-40	Profox software <i>Download cable sold separately.</i>
1053371	Download cable
1019441	nVision software <i>Download cable sold separately.</i>
100-0433-00	Battery door

See the replacement parts catalog for additional replacement parts.



### 920 and 920M hand-held oximeters

The 920 and 920M oximeters are small, easy-to-use handheld devices that deliver on-the-spot oximetry readings; perfect for home use. The 920M has the added capability of storing up to 18 hours of trend memory. Both oximeters are intended for spot check use and can be connected to an external printer.



# Oximeters

Part number	
920MP	920M Plus hand-held oximeter with memory <i>Includes 936 adult finger clip sensor, four AA batteries, and a manual.</i>
1028009	2500A hand-held oximeter downloadable with alarms

Accessories	
Software	
H920M-40	Profox software <i>Download cable sold separately.</i>
1019441	nVision software <i>Download cable sold separately.</i>
Miscellaneous	
920MPC-NA	920MP rechargeable battery system, wall mount <i>Includes 920MPC charger stand, NiMH rechargeable battery, power supply, and an operator manual.</i>
1053371	Download cable
920MPC	Charger stand only, replacement
1007355	920MP NiMH rechargeable battery
1007332	Carrying case

See the replacement parts catalog for additional replacement parts.



## 920M Plus hand-held oximeter with memory

The 920M Plus handheld pulse oximeter provides on-the-spot oximetry readings, and is designed to be easy to use for both healthcare professionals and patients. The 920M Plus features overnight recordings (72 hours of memory), is compatible with all Respironics sensors, uses disposable or rechargeable batteries, and interfaces with Profox and nVISION software.



## 2500A hand-held oximeter with memory and alarms

The 2500 handheld oximeter is a small, versatile digital pulse oximeter for accurate assessment of blood oxygen saturation and pulse rate. It is ideal for ambulatory monitoring and overnight recordings. The 2500A has two-button operation, is compact, has three power options, and features visual and audible alarms. Compatible with nVision and Profox software.

# Oximeters

## Part number

1062232	GO <sub>2</sub> finger oximeter
---------	---------------------------------

## Accessories

1062965	Carrying case
1062966	Lanyard



### GO<sub>2</sub> pulse oximeter

The quality construction and high accuracy of the GO<sub>2</sub> finger oximeter allows your patients to take greater control of their lives while creating a new opportunity for your business. Oxygen patients now have a quality alternative in an accurate, affordable, and portable oximeter allowing them the freedom they need to go about their activities. You have the opportunity to offer your patients a personal oximeter that provides accuracy comparable to professional oximeters. Your patients will benefit from GO<sub>2</sub>'s reliable and accurate oxygen saturation and heart rate monitoring information, which allows them to manage both their oxygen and their activity level.

# Oximeters

## Part number

1078606	WristOx <sub>2</sub> 3150, pulse oximeter with standard accessories (reusable soft sensor, wrist bands, batteries, CD)
---------	--

## Accessories

1078747	Medium soft sensor, 8 pin (12"/0.3 m)
1078749	Adult flex sensor, 8 pin (12"/0.3 m)
1079107	Adult finger clip sensor, 8 pin (12"/0.3 m)
1084509	Sensor adapter, WristOx <sub>2</sub> , 3150 (8 to 9 pin)
1078748	Sensor adapter, WristOx 3100 (9 to 8 pin)
1078750	USB PC download cable for WristOx2
1079118	Wrist band, reusable (10")
1091490	WristOx <sub>2</sub> , wrist band (8")
1091489	WristOx <sub>2</sub> , wristband (6")



### WristOx<sub>2</sub> ambulatory digital pulse oximeter

WristOx<sub>2</sub> is a wrist-worn pulse oximeter that has been developed to deliver highly accurate readings and help reduce the difficulties that can be associated with overnight oximetry studies. Key design features include an enhanced memory, the use of standard batteries, the addition of wireless data transfer, and more.

# Oxygen accessories

Part number		HCPCS
<b>Carrying cases</b>		
900-107-02	Black carrying case <i>For use with the EverGo oxygen concentrator.</i>	
900-250	Accessory case (pictured) <i>For use with the EverGo oxygen concentrator.</i>	
1018019	GoLox carrying case	
1062965	GO <sub>2</sub> carrying case	
1019439	WristOx ambulatory digital pulse oximeter carrying case	
1007332	Oximeter carrying case <i>For use with 920, 920M, and 920MP oximeters.</i>	
<b>Cannulas</b>		
1025655	Cannula, 4' <i>For use with ePOD.</i>	A4615*
1014653	High flow starter kit <i>Includes cannula, 25' tubing, connector. For use with M10, M10600, and M10605 oxygen concentrators.</i>	4615**, A4616**
1014654	High flow nasal cannula <i>For use with M10, M10600, and M10605 oxygen concentrators.</i>	A4615**
1017768	High flow oxygen tubing, 25' <i>For use with M10, M10600, and M10605 oxygen concentrators.</i>	A4616**
1014655	Oxygen tubing connector <i>For use with M10, M10600, and M10605 oxygen concentrators.</i>	


\*Not separately Medicare reimbursable.

\*\*Accessories are not separately billable to Medicare, except when used with patient-owned equipment. Other payers may reimburse separately for these items. Check with the specific payer.

Part number	
<b>Extension cables</b>	
927	Extension cable, 3' <i>For use with 920, 920M, and 930 oximeters.</i>
1008130	Extension cable, 10' (pictured) <i>For use with 920 and 920M oximeters.</i>
<b>Filters</b>	
1038831	Compressor intake filter <i>For use with EverFlo oxygen concentrators.</i>
900-150	Air inlet filters, 6 per pk. (pictured) <i>For use with EverGo oxygen concentrators.</i>
520	Micro disk filter <i>For use with EverFlo oxygen concentrators.</i>
1017796	Inlet filter replacement kit <i>For use with M10, M10600, and M10605 oxygen concentrators.</i>
1015581	Kit, grommet, 10 pk.
H618	Inlet pre-filter replacement kits, 6 per pk. <i>For use with all Millennium oxygen concentrators.</i>
H619	Air inlet filter replacement kits, 6 per pk. <i>For use with M600 and M605 oxygen concentrators.</i>
H620	Inlet filter replacement kits, 6 per pk. <i>For use with 600, 605, M600, and M605 oxygen concentrators.</i>
H621	Microdisk filter replacement kits, 6 per pk. <i>For use with all Millennium oxygen concentrators.</i>
H652	Filter cover for compressor kit <i>For use with 600, 605, M600, and M605 oxygen concentrators.</i>
H653	Filter felt for compressor kits, 6 per pk. <i>For use with 600, 605, M600, and M605 oxygen concentrators.</i>



1008130



900-150



528

# Oxygen accessories

Part number	
Flow meters	
H537-1	Stand-alone low flow meter kit <i>For use with EverFlo oxygen concentrators.</i>
H644	Cabinet integrated pediatric flow meter kit <i>For use with 600, 605, and M600 (non-OPI units only) oxygen concentrators.</i>
528	Locking flow meter kit (pictured) <i>For use with 600, 605, M600, and M605 oxygen concentrators.</i>
1082784	5 LPM flow meter replacement kit <i>For use with 600, 605, M600, and M605 oxygen concentrators.</i>
Gaskets	
1007715	Gasket replacement kit <i>For use with all Millenniums.</i>
1007716	Gasket replacement tool kit <i>For use with 600, 605, M600, and M605 oxygen concentrators.</i>

Part number	
Filter caps	
H649	Caster brake kit <i>For use with EverFlo oxygen concentrators.</i>
900-100	Mobile cart (pictured) <i>For use with the EverGo and SimplyGo oxygen concentrators.</i>
900-102	Rechargeable li ion battery (pictured) <i>For use with the EverGo oxygen concentrator.</i>
900-103	Battery charger/recalibrator (pictured) <i>For use with the EverGo oxygen concentrator.</i>
900-121	Airline connector <i>For use with the EverGo oxygen concentrator.</i>
920MPC	Charger stand <i>Power cord not included. For use with 920M.</i>
100-0339-00	External spring, single or 30 per case, (pictured) <i>For use with 950 oximeter.</i>
100-0338-00	Gripping insert <i>For use with 950 oximeter.</i>
1007355	920MP NiMH rechargeable battery
1017984	Shoulder strap <i>For use with GoLox.</i>
1062966	Lanyard, black, 20" <i>For use with GO2 oximeter</i>
H654	Air inlet retainer kits, 10 per pk. <i>For use with 600 and 605 oxygen concentrators.</i>
1001684	Flow signal conditioner kit
1019440	Reusable wristbands, 10 per pk. <i>For use with WristOx.</i>
1019438	Reusable wristband extenders, 5", 10 per pk. <i>For use with WristOx.</i>
1015581	Silicone grommets, 10 per pk. <i>For use with M10600, and M10605 oxygen concentrators.</i>



900-100



900-102




900-103



# Oxygen accessories

Part number	
<b>Oxygen attachments</b>	
7050	Tee adapter, oxygen
7314	Tubing adapters, 22mm x 22mm, 6 per pk.
300065	Tubing coupler
312010	O <sub>2</sub> enrichment attachments, 10 per pk.
312710	O <sub>2</sub> enrichment attachment
1039642	Humidifier connector tube
<b>Oxygen valves</b>	
302418	Pressure valve
302423	Pressure valves, 5 per pk.
1000348	Oxygen valve
<b>Power cords</b>	
900-104	AC power cord <i>For use with the EverGo oxygen concentrator.</i>
900-101	DC power cord <i>For use with the EverGo (serial number &lt; 109500)</i>
900-120	DC power cord (pictured) <i>For use with the EverGo (serial number &gt;= 109500)</i>
H617	Power cord replacement kit, 60Hz <i>For use with all Millenniums.</i>
<b>Power supplies and battery packs</b>	
900-105	AC power supply (pictured) <i>For use with the EverGo oxygen concentrator.</i>
<b>Sensors: reusable sensors</b>	
1037296	Soft sensor, small, 1M
1037297	Soft sensor, medium, 1M
1037298	Soft sensor, large, 1M
1040731	Soft sensor, medium <i>For use with WristOx.</i>

Part number	
<b>Sensors: disposable sensors</b>	
1037302	Disposable sensors, pediatric small, 24 per pk.
1037300	Disposable sensors, pediatric, 24 per pk.
1037301	Disposable sensors, infant, 24 per pk.
1037299	Disposable sensors, adult, 24 per pk.
953	Adult sensor, 44 lbs. or more <i>Includes Flex sensor, 38mm, straight cable, 25 Hydrogel tape strips, and Posey wraps, 1 meter cable.</i>
953-11	Adult sensor, 3 meter cable, 8'
934	Infant sensor, large toe, 4 to 44 lbs. <i>Includes Flex sensor, 34mm, with FlexiWrap tape strips.</i>
935	Pediatric small sensor, outstep of foot, 4 lbs. or less <i>Includes Flex Sensor, 38mm, right angle cable, with 25 FlexiWrap tape strips and 12 Posey wraps.</i>
936	Adult finger clip sensor, index finger, 66 lbs. or more
936P	Pediatric finger clip sensor, index finger, 18 lbs. or more
948	Adult ear clip sensor, lobe of ear, 88 lbs. or more
951	Reflectance sensor middle forehead, 66 lbs. or more
1019522	Finger clip sensor, adult <i>For use with WristOx.</i>
<b>Sensor cables</b>	
1019444	Sensor adapter cable <i>For use with WristOx.</i>
938	Patient cable; connects sensor assembly to oximeter (pictured) <i>For use with 930 Oximeter.</i>
<b>Sensor tape</b>	
940	Double-sided tape for 934 sensor, 100 per pk.
941	Double sided tape for 935 and 953, 100 per pk.
942	Overwrap tape strips, 100 per pk.
947	Posey wrap for 935 and 953 sensor, 12 per pk. (pictured)
954A	Adult FlexiWrap tape strips for 953 sensor, 25 per pk.
954I	Infant FlexiWrap tape strips for 934 sensor, 25 per pk.
954N	Pediatric small FlexiWrap tape strips for 935 sensor, 25 per pk.



100-0339-00



900-120



900-105



938



947



# Masks

Our complete line of nasal, full-face, minimal contact, and pediatric masks reflect our commitment to innovation and meeting patient and business needs.

Not all products are available in all countries. Check with Philips Respironics customer service or your local distributor.



# Nuance

A new standard. A new choice. Nuance gel pillows platform.

We're setting a new standard in nasal pillows masks. Philips Respironics, the world leader in gel cushion technology, is bringing the comfort and sealing power of gel to nasal pillows. As the latest in our groundbreaking new line of masks, Nuance gel pillows represents the new choice for better sleep therapy for clinicians, home care providers and patients.

As the first and only gel pillows platform, Nuance is designed to help patients succeed by providing a more comfortable and familiar feel. With the greater comfort and seal of gel pillows, the ability to fit many different size nostrils and choice of a fabric or gel frame, Nuance is truly a new standard in nasal pillows masks.

#### Key features

- Choice between fabric or gel frame for comfort with fewer red marks
- Non-slip headgear stays in place throughout the night
- Gel pillows provides greater comfort and seal
- Lightweight, flexible tubing
- Reduces nasal irritation
- Conforms to different size nostrils



Part number		HCPGS
<b>Nuance</b>		
1105160	Nuance gel pillows mask, with headgear	A7034,A7035
1106194	Nuance gel pillows mask, without headgear	A7034
<b>Nuance Pro</b>		
1105167	Nuance Pro gel pillows mask, with headgear	A7034,A7035
1106195	Nuance Pro gel pillows mask, without headgear	A7034
<b>Replacement parts</b>		
1105176	Nuance headgear	A7035
1105405	Nuance RS fabric frame	
1105177	Nuance fabric frame	
1105178	Nuance Pro headgear	A7035
1105179	Nuance Pro frame	
1105173	Nuance/Nuance Pro small gel cushion	A7033
1105174	Nuance/Nuance Pro medium gel cushion	A7033
1105175	Nuance/Nuance Pro large gel cushion	A7033
1105180	Nuance/Nuance Pro swivel tube with exhalation	A7045



# Nuance. A new choice.

Designed around patients.

## Gel frame

- Gel-padded frame and non-slip headgear holds mask in place
- Reduces need for re-adjustment
- High-performance look



## Fabric frame

- All fabric frame and headgear enhances patient comfort
- Soft fabric experience for better sleep
- Fewer red marks



To learn more, visit

<http://www.theanswer.philips.com/nuance/>



# Wisp

Delivering a combination of comfort, performance and style.

Developing solutions that meet the needs of a rapidly changing industry calls for fresh thinking and collaboration. At Philips Respironics we traveled around the world asking clinicians, homecare providers, and patients how we could improve their experience.

Clinicians wanted us to resolve the tradeoff between nasal and pillows masks, homecare providers wanted a solution that could reduce operational complexity by getting more patients served from one package, and patients wanted a better overall mask experience that improved their lifestyle.

When you put it all together, the answer is Wisp.

## Key features

- Compact design for a natural fit
- Soft and comfortable for a restful night's sleep
- Open field of vision for a feeling of independence
- Quick and easy to adjust for a superior seal and personalized fit
- Contemporary frame options to complement your style



Part number		HCPCS
<b>Wisp, fabric frame</b>		
1094051	Wisp mask with fabric frame, with headgear	A7034,A7035
1094056	Wisp mask with fabric frame, without headgear	A7034
<b>Wisp, clear frame</b>		
1094050	Wisp mask with clear frame, with headgear	A7034,A7035
1094055	Wisp mask with clear frame, without headgear	A7034
<b>Wisp replacement parts</b>		
1094086	Petite cushion	A7033
1094087	Small/medium cushion	A7033
1094088	Large cushion	A7033
1094080	Clear frame	
1094081	Fabric frame	
1094082	Wisp headgear	A7035
1094078	Wisp headgear, large size	A7035
1094083	Wisp headgear clips	



All three sizes in one package: P, SM, L




Compact design comfortably fits more than **98%** of patients



Patients said they had a **better quality of sleep** with Wisp over the leading pillows and nasal masks



Open field of vision for a feeling of **greater independence** – patients can read and even wear their glasses



Part number		HCPCS
<b>TrueBlue with headgear</b>		
1071800	Petite	A7034,A7035
1071801	Small	A7034,A7035
1071803	Medium	A7034,A7035
1071804	Medium wide	A7034,A7035
1071805	Large	A7034,A7035
<b>TrueBlue without headgear</b>		
1071806	Petite	A7034
1071807	Small	A7034
1071808	Medium	A7034
1071809	Medium wide	A7034
1071810	Large	A7034
<b>TrueBlue with headgear DuoPack (includes two same-size cushion sets)</b>		
1071818	Petite	A7034,A7035,A7032
1071819	Small	A7034,A7035,A7032
1071820	Medium	A7034,A7035,A7032
1071821	Medium wide	A7034,A7035,A7032
1071822	Large	A7034,A7035,A7032
<b>TrueBlue with headgear FitPack (includes cushions in two sizes)</b>		
1071816	Small and medium	A7034,A7035
1071817	Medium wide and large	A7034,A7035
<b>TrueBlue replacement parts</b>		
1071861	Petite cushion and flap	A7032
1071862	Small cushion and flap	A7032
1071863	Medium cushion and flap	A7032
1071864	Medium wide cushion and flap	A7032
1071865	Large cushion and flap	A7032
1071875	TrueBlue headgear	A7035
1071874	TrueBlue headgear, reduced size	A7035
1071888	TrueBlue gel forehead pad	
1071877	TrueBlue talon clips (pack of 2)	



**TrueBlue**

TrueBlue brings our best technologies together in one mask. Designed to push the limits of performance and fit, this mask sets the gold standard for a good night's sleep. The TrueBlue gel nasal mask with Auto Seal technology is designed to deliver a higher degree of comfort, stability, and freedom of movement with minimal adjustments. This intuitive gel mask is designed for quicker set-up and fitting, fewer callbacks, and easier lab titration.

## ComfortGel Blue

Part number		HCPCS
<b>ComfortGel Blue with headgear</b>		
1070040	Petite	A7034,A7035
1070039	Small	A7034,A7035
1070038	Medium	A7034,A7035
1070037	Large	A7034,A7035
<b>ComfortGel Blue without headgear</b>		
1070044	Petite	A7034
1070043	Small	A7034
1070042	Medium	A7034
1070041	Large	A7034
<b>ComfortGel Blue with headgear DuoPack (includes two same-size cushion sets)</b>		
1070053	Small and medium	A7034,A7035,A7032
<b>ComfortGel Blue replacement parts</b>		
1070108	Petite cushion and flap	A7032
1070107	Small cushion and flap	A7032
1070106	Medium cushion and flap	A7032
1070105	Large cushion and flap	A7032
1070104	Petite cushion	A7032
1070103	Small cushion	A7032
1070102	Medium cushion	A7032
1070101	Large cushion	A7032
1033678	Premium headgear	A7035
1059236	RS Premium headgear (smaller)	A7035
1040114	Forehead pad	
1009723	ComfortGel Blue Quick Clips (2 per pack)	



### ComfortGel Blue

Staying true to our acclaimed blue gel masks, we've retained the globally popular comfort and flexibility features you expect, while enhancing patient comfort even further. With ComfortGel Blue, we've minimized pressure points, reduced noise, and redirected exhalation air flow away from the bed partner. These and other enhancements may help to create a better patient experience and encourage greater compliance.

Part number		HCPCS
1050000	EasyLife nasal mask with headgear, petite	A7034,A7035
1050001	EasyLife nasal mask with headgear, small	A7034,A7035
1050002	EasyLife nasal mask with headgear, medium	A7034,A7035
1050003	EasyLife nasal mask with headgear, medium wide	A7034,A7035
1050004	EasyLife nasal mask with headgear, large	A7034,A7035
1050005	EasyLife DuoPack with headgear, petite (two cushion sets*)	A7034,A7035,A7032
1050078	EasyLife FitPack with headgear, S, M cushion sets*	A7034,A7035
1050079	EasyLife FitPack with headgear, MW, L cushion sets*	A7034,A7035
1050006	EasyLife DuoPack with headgear, small (two cushion sets*)	A7034,A7035,A7032
1050007	EasyLife DuoPack with headgear, medium (two cushion sets*)	A7034,A7035,A7032
1050008	EasyLife DuoPack with headgear, medium wide (two cushion sets*)	A7034,A7035,A7032
1050009	EasyLife DuoPack with headgear, large (two cushion sets*)	A7034,A7035,A7032

Accessories		HCPCS
1050005	EasyLife without headgear, petite	A7034
1050090	EasyLife cushion set*, petite	A7032
1050011	EasyLife without headgear, small	A7034
1050012	EasyLife without headgear, medium	A7034
1050013	EasyLife without headgear, medium wide	A7034
1050014	EasyLife without headgear, large	A7034
1050091	EasyLife cushion set*, small	A7032
1050092	EasyLife cushion set*, medium	A7032
1050093	EasyLife cushion set*, medium wide	A7032
1050094	EasyLife cushion set*, large	A7032
1050095	EasyLife headgear	A7035
1050096	EasyLife sizing gauge	

\*A cushion set includes one inner seal cushion and one outer support cushion.

Other headgear can be used with this mask including premium headgear, RS premium headgear, deluxe headgear, and small deluxe headgear.



**EasyLife nasal mask**

With its patented Auto Seal technology, the EasyLife nasal mask nearly eliminates the need for manual mask adjustments. Its innovative, lightweight design features a unique dual-cushion construction - the inner cushion creates an instant, self-adjusting seal, while the outer cushion provides comfortable support. The only manual steps necessary are intuitive headgear adjustments, making fitting as easy as possible.

Part number		HCPCS
1007963	ComfortClassic nasal mask with deluxe headgear, small	A7034,A7035
1007964	ComfortClassic nasal mask with deluxe headgear, medium	A7034,A7035

Accessories		HCPCS
1007966	ComfortClassic without headgear, small	A7034
1007965	ComfortClassic without headgear, medium	A7034
152033	Port caps, 10 per pk.	
1002464	ComfortClassic gel forehead pad replacement kits, 2 per pk.	
1008651	Foam forehead spacers, 5 per pk.	
1004395	Elbow with exhalation	A7045
1033678	Premium headgear with EZPeel tabs	A7035
1059236	RS premium headgear	A7035
1002800	Deluxe headgear	A7035
1018755	Small deluxe headgear	A7035
1002720	Split washer kits, 5 per pk.	
1012041	Sizing gauge	

*Other headgear can be used with this mask including premium headgear, RS premium headgear, deluxe headgear, and small deluxe headgear.*

*Accessories are not separately billable to the Medicare program, with the exception of replacement circuits. Consult your local Medicare DME MAC local coverage determination (LCD) for details. Commercial insurance plans may provide separate coverage for accessories. Consult insurance company policy manual for more information.*



**ComfortClassic nasal mask**

ComfortClassic features a dual-layer silicone cushion with an outer membrane that fills with air to create a comfortable seal. The thicker inner section provides stability. A dual-position gel spacer molds to the patient’s forehead for even greater comfort and stability, and a “silent” exhalation port gives you a mask that may help to improve patient compliance.



## Profile Lite

Part number		HCPCS
1004086	Profile Lite mask with deluxe headgear, petite	A7034,A7035
1004087	Profile Lite mask with deluxe headgear, small	A7034,A7035
1004088	Profile Lite mask with deluxe headgear, medium small	A7034,A7035
1004089	Profile Lite mask with deluxe headgear, medium	A7034,A7035
1004110	Profile Lite mask with deluxe headgear, medium wide	A7034,A7035
1004111	Profile Lite mask with deluxe headgear, large	A7034,A7035
1004112	Profile Lite mask with deluxe headgear, large narrow	A7034,A7035
1047339	Profile Lite youth mask with deluxe headgear	A7034,A7035

All units include an exhalation port and appropriately sized gel spacer.

Accessories		HCPCS
1002849	Profile Lite without headgear, petite	A7034
1002850	Profile Lite without headgear, small	A7034
1002370	Profile Lite without headgear, medium small	A7034
1002371	Profile Lite without headgear, medium	A7034
1002372	Profile Lite without headgear, medium wide	A7034
1002373	Profile Lite without headgear, large	A7034
1002374	Profile Lite without headgear, large narrow	A7034
1004395	Elbow with exhalation	A7045
152033	Port caps, 10 per pk.	
1004292	Small gel retrofit forehead spacer	
1004291	Large gel retrofit forehead spacer	
1008869	Extra large gel retrofit forehead spacer	
1009857	Gel spacer retrofit assortment pack <i>Includes 3 small, 1 large, 3 extra large</i>	
1002800	Deluxe headgear	A7035
1018755	Small deluxe headgear	A7035
302429	Sizing gauge	

Other headgear can be used with this mask including premium headgear, RS premium headgear, deluxe headgear, and small deluxe headgear.



### Respironics Profile Lite nasal gel mask and Respironics Profile Lite youth mask

The soft, light cushion is extremely comfortable right out of the box. The gel cushion and forehead pad gently rest against the patient's face for a soft, pillowy fit. The Profile Lite is fully customizable. For a more personalized seal, a 10-minute customization process molds the mask to the individual's facial contours.



# Amara

Performance, appeal, and simplicity at its best

With gel and silicone options and our one “click” design, Amara delivers performance and ease of use for your patients. With the addition of a gel cushion option, Philips Respironics has made the ingeniously-simple Amara full face mask even more brilliant. Your patients can now choose the type of comfort they prefer in their cushion with gel or silicone options.

Amara offers an excellent seal with comfort and ease-of-use advantages to current traditional leading full face offerings. Both the gel and silicone options are smaller and lighter weight than the leading full face mask, and have 60% fewer parts. Amara is lightweight with a simple one “click” design and cushion options that fits over 95% of faces.

## Key features

- Easy to fit, assemble and use
- Smallest/lightest leading traditional full-face mask
- 60% fewer parts than the leading full-face mask
- Proven to fit 95% of patients' faces

To learn more, visit

<http://www.theanswer.philips.com/amara/>

**One “click” cushion assembly**

A cushion design that reduces pressure points.

**Quick-fit headgear clips**

Reduces the need for continual headgear adjustment.

**Fine-glide forehead adjuster**

Ten selectable positions for fine-tuning seal and fit.

**Lightweight minimal footprint frame**


Smaller and lighter than any other leading traditional full-face mask on the market.

**360° rotating elbow**

Quiet-exhalation micro ports with quick release swivel.

**Four-part modular design**

Modular design for quick assembly and cleaning, with interchangeable cushions available in four sizes.



Silicone		
Part number		HCPCS
<b>Amara with headgear</b>		
1090200	Petite	A7030,A7035
1090201	Small (reduced size frame and headgear)	A7030,A7035
1090226	Medium (reduced size frame)	A7030,A7035
1090227	Large (reduced size frame)	A7030,A7035
<b>Amara without headgear</b>		
1090210	Petite	A7030
1090211	Small (reduced size frame)	A7030
1090215	Medium (reduced size frame)	A7030
1090216	Large (reduced size frame)	A7030
<b>Amara replacement parts</b>		
1090290	Petite cushion	A7031
1090292	Small cushion	A7031
1090293	Medium cushion	A7031
1090294	Large cushion	A7031
1090288	Frame, reduced size (accepts all cushions)	
1090289	Frame, standard (accepts all cushions)	
1090296	Amara headgear, reduced size	A7035
1090297	Amara headgear, standard	A7035
1040114	Forehead pad	
1071877	Talon headgear clips	
1081899	Entrainment valve swivel	



Available in four sizes: P, S, M, L

Amara starter kits (includes gel mask with same size silicone cushion)		
1106618	Petite	A7030,A7035,A7031
1106619	Small	A7030,A7035,A7031
1106620	Medium	A7030,A7035,A7031
1106621	Large	A7030,A7035,A7031

Gel		
Part number		HCPCS
<b>Amara with headgear</b>		
1090400	Petite (reduced size frame and headgear)	A7030,A7035
1090401	Small (reduced size frame and headgear)	A7030,A7035
1090405	Medium (reduced size frame)	A7030,A7035
1090406	Large (reduced size frame)	A7030,A7035
<b>Amara without headgear</b>		
1090410	Petite (reduced size frame)	A7030
1090411	Small (reduced size frame)	A7030
1090415	Medium (reduced size frame)	A7030
1090416	Large (reduced size frame)	A7030
<b>Amara replacement parts</b>		
1090490	Petite cushion	A7031
1090492	Small cushion	A7031
1090493	Medium cushion	A7031
1090494	Large cushion	A7031
1090296	Amara headgear, reduced size	A7035
1090297	Amara headgear, standard	A7035
1040114	Forehead pad	
1071877	Talon headgear clips	
1081899	Entrainment valve swivel	



Available in four sizes: P, S, M, L

## ComfortGel Blue Full

Part number		HCPCS
1040135	ComfortGel Blue Full mask with premium headgear, small	A7030, A7035
1040136	ComfortGel Blue Full mask with premium headgear, medium	A7030, A7035
1040137	ComfortGel Blue Full mask with premium headgear, large	A7030, A7035
1062004*	ComfortGel Blue Full mask with premium headgear, extra large	A7030, A7035
1002854	ComfortGel Blue Full mask DuoPack with premium headgear, small plus small cushion	A7030, A7035, A7031
1002855	ComfortGel Blue Full mask DuoPack with premium headgear, medium plus medium cushion	A7030, A7035, A7031
1002856	ComfortGel Blue Full mask DuoPack with premium headgear, large plus large cushion	A7030, A7035, A7031
1062005*	ComfortGel Blue Full mask DuoPack with premium headgear, extra large plus extra large cushion	A7030, A7035, A7031

Accessories		HCPCS
1040568	ComfortGel Blue Full mask without headgear, small	A7030
1040569	ComfortGel Blue Full mask without headgear, medium	A7030
1040570	ComfortGel Blue Full mask without headgear, large	A7030
1062003*	ComfortGel Blue Full mask without headgear, extra large	A7030
1051032	ComfortGel Blue Full mask without headgear, DuoPack, small plus small cushion	A7030, A7031
1051033	ComfortGel Blue Full mask without headgear, DuoPack, medium plus medium cushion	A7030, A7031
1051034	ComfortGel Blue Full mask without headgear, DuoPack, large plus large cushion	A7030, A7031
1040798	ComfortGel Blue Full cushion, small	A7031
1040799	ComfortGel Blue Full cushion, medium	A7031
1040800	ComfortGel Blue Full cushion, large	A7031
1062030*	ComfortGel Blue Full cushion, extra large	A7031
1040795	SST flap only, small	A9999
1040796	SST flap only, medium	A9999
1040797	SST flap only, large	A9999
1062031*	SST flap only, extra large	A9999
1044861	ComfortGel Blue Full cushion and SST flap, small	A7031
1044862	ComfortGel Blue Full cushion and SST flap, medium	A7031
1044863	ComfortGel Blue Full cushion and SST flap, large	A7031
1062032*	ComfortGel Blue Full cushion and SST flap, extra large	A7031
1040138	ComfortGel Blue Full premium headgear	A7035
1040114	Premium forehead pad	
1049398	Entrainment valve swivel	
1040139	Ball-and-socket Quick Clips, blue, 2 per pk.	
1010872	Port caps, 5 per pk.	
1049507	Sizing gauge	



### ComfortGel Blue Full full-face mask

Combines an all-new design with some very familiar – and popular – components. Our comfortable forehead cushion and enhanced StabilitySelector combine with our popular blue gel replaceable cushion to create a level of comfort and stability no other full-face mask can approach. The unique sure seal technology (SST) provides a membrane-thin layer over the gel cushion to follow the contours of the face better, allowing wearers to quickly achieve the comfortable seal they need.

*Other headgear can be used with this mask including premium headgear, RS premium headgear, deluxe headgear, and small deluxe headgear.*



Part number		HCPCS
<b>FitLife with headgear</b>		
1060801	Small	A7030,A7035
1060802	Large	A7030,A7035
1089991	Extra large	A7030,A7035
<b>FitLife without headgear</b>		
1072530	Small	A7030
1072531	Large	A7030
1089992	Extra large	A7030
<b>FitLife replacement parts</b>		
1060807	FitLife headgear, small	A7035
1060808	FitLife headgear, large	A7035
1089990	FitLife headgear, extra large	A7035
1060913	FitLife snap clips, 2 per pk.	
7041	Mask swivel	
1060809	Entrainment elbow	A7045



### FitLife

A solution for patients who have difficulty wearing conventional masks or who simply prefer a clear line of sight. FitLife addresses a number of patient issues. Facial hair, dentures, nasal bridge breakdown, facial irregularities, and claustrophobia are issues that can greatly impact the ability to establish effective therapy and comfort. FitLife quickly establishes a robust seal around the less pressure-sensitive perimeter of the face, while maintaining an unrestricted view.

Part number		HCPCS
<b>GoLife for Women</b>		
1082577	FitPack with headgear Includes petite, small, medium cushions	A7034,A7035
1082579	FitPack without headgear Includes petite, small, medium cushions	A7034
<b>GoLife for Women subassembly</b>		
1082578	Without headgear, without cushions	
1082583	Subassembly only, no swivel tubing	
<b>GoLife for Men</b>		
1073114	FitPack with headgear Includes small, medium, large cushions	A7034,A7035
1073329	FitPack without headgear Includes small, medium, large cushions	A7034
1093275	FitPack without headgear (small frame) Includes small, medium, large cushions	A7034
<b>GoLife for Men subassembly</b>		
1073328	Without headgear, without cushions	
1093274	Small frame, without headgear, without cushions	
1078408	Subassembly only, no swivel tubing	
1093229	Small frame, subassembly only, no swivel tubing	
<b>GoLife replacement parts</b>		
1073332	Petite cushion	A7033
1073333	Small cushion	A7033
1073334	Medium cushion	A7033
1073335	Large cushion	A7033
1082581	GoLife for Women headgear	A7035
1073336	GoLife for Men headgear	A7035
1073337	Tube management clip, 2 per pk.	
1078407	Swivel tubing (18"/45.7cm.) with exhalation port	A7045



**GoLife**

GoLife for Men is designed specifically to fit the larger contours of the male facial structure, while GoLife for Women is developed with a distinctly smaller frame and headgear to more closely fit feminine facial dimensions. Both feature straightforward headgear and self-adjusting, optimally angled nasal pillows for a one-step fit. Quite simply, these unique GoLife masks help to take the complexity out of achieving successful nasal pillows therapy for both male and female patients.

Part number		HCPCS
1036800	OptiLife FitPack with headgear and chin support band size P, S, M, L pillows cushions	A7034, A7035
1036801	OptiLife FitPack/DuoPack with headgear and chin support band size P, S, M pillows cushions and S, M CradleCushions	A7034, A7035, A7032
1036802	OptiLife FitPack/DuoPack with headgear and chin support band size M, L pillows cushions and M, L, LN CradleCushions	A7034, A7035, A7032
1036818	OptiLife mask sub-assembly with headgear and chin support band <i>No cushions</i>	A7035*

Accessories		HCPCS
1036819	OptiLife sub-assembly <i>No cushions or headgear</i>	*
1036804	OptiLife FitPack, mask only size P, S, M, L pillows cushions	A7034
1036838	OptiLife pillows cushion, petite	A7033
1036839	OptiLife pillows cushion, small	A7033
1036840	OptiLife pillows cushion, medium	A7033
1036841	OptiLife pillows cushion, large	A7033
1036846	OptiLife CradleCushion, small	A7032
1036847	OptiLife CradleCushion, medium	A7032
1036848	OptiLife CradleCushion, large narrow	A7032
1036849	OptiLife CradleCushion, large	A7032
1036850	Headgear <i>Includes chin support band</i>	A7035
1039967	Tube management clips, 2 per pk.	
1032007	Respironics pillows sizing gauges, 5 per pk.	
1031154	Respironics pillows sizing gauges, 20 per pk.	

\*Requires addition of mask cushion for HCPCS A7034.



**OptiLife nasal mask**

Features a headgear design with an integrated chin support band that allows it to be put on and adjusted with little effort – and no buckles.



CradleCushion



Pillows cushion

## Accessories/U.S. HCPCS codes

Part number	
<b>Accessories</b>	
302175	Chin strap
1012911	Premium chin strap
302425	Deluxe chin strap
1032907	Lightweight CPAP tubing (6 ft.)
7212	Standard CPAP tubing (10 ft.)
1036425	Sizing gauge ring
<b>Codes for nasal masks</b>	
A7034	Nasal mask only
A7035	Headgear
A7032	Nasal cushion replacement
<b>Codes for full-face masks</b>	
A7030	Full-face mask only
A7035	Headgear
A7031	Full-face cushion replacement
<b>Codes for nasal pillows masks</b>	
A7034	Nasal mask only
A7035	Headgear
A7033	Nasal pillow replacement
A7032	Nasal cushion replacement (OptiLife CradleCushion)
A7045	Exhalation port with swivel (GoLife PN 1078407)



# Sleep diagnostics

Whether you're considering new PSG systems, home sleep testing systems, titration systems, sensors or other diagnostic tools, we can help you with our full line of in-lab and portable diagnostic products, services and solutions. We offer unlimited software licensing, core software updates, starter sensor kits, and more to help you realize long-term benefits through operational efficiencies.

Not all products are available in all countries. Check with Philips Respironics customer service or your local distributor.



# Alice 6 LDx sleep systems

Pair the LDx base station with either the LDxS or LDxN headbox for a full-featured, easy-to-understand system that allows staff to focus more on patients and less on equipment. The LDxS is the perfect front-line headbox for most sleep center needs. Use the LDxN when you need to increase from 19 EEG inputs on the LDxS to a full 32 EEG inputs. The LDx base station works with either headbox, so it's easy to have the right capability where, and when, you need it.



**Alice 6 LDx base system,  
North America**  
**1063314**

Includes:

- Alice 6 LDx base station
- AC power cord
- Auto-switching power supply
- Alice 6 communication cable
- Alice 6 user guide
- Headbox shoulder strap (headbox not included)
- Unlimited site-licensing
- 2-year warranty
- Software upgrades
- Starter zRIP DuraBelt kit



**LDxS Headbox**  
**1063552** North America  
**1063353** International



**Alice 6 LDx base system,  
International**  
**1063315**

Includes:

- Alice 6 LDx base station
- AC power cord
- Auto-switching power supply
- Alice 6 communication cable
- Alice 6 user guide
- Headbox shoulder strap (headbox not included)
- Unlimited site-licensing
- 2-year warranty
- Software upgrades
- Starter zRIP DuraBelt kit



**LDxN Headbox**  
**1063476** North America  
**1063547** International

To learn more visit [www.philips.com/sleepdx](http://www.philips.com/sleepdx).

# Alice 6 LDE sleep system

The Alice 6 LDE diagnostic system provides your busy lab with the base set of channels needed to meet AASM standards while allowing staff to focus more on patients and less on equipment.



**Alice 6 LDE base system,  
North America**  
**1063311**

Includes:

- Alice 6 LDE base station
- AC power cord
- Auto-switching power supply
- Alice 6 communication cable
- Alice 6 user guide
- Headbox shoulder strap (headbox not included)
- Unlimited site-licensing
- 1-year warranty
- Software upgrades
- Starter zRIP DuraBelt kit



**Headbox, North America**  
**1063549**

Alice 6 LDE headbox,  
North America



**Alice 6 LDE base system,  
International**  
**1063312**

Includes:

- Alice 6 LDE base station
- AC power cord
- Auto-switching power supply
- Alice 6 communication cable
- Alice 6 user guide
- Headbox shoulder strap (headbox not included)
- Unlimited site-licensing
- 1-year warranty
- Software upgrades
- Starter zRIP DuraBelt kit



**Headbox, International**  
**1063550**

Alice 6 LDE headbox,  
International

# Computers/monitors



1013805



1043655



1043167



1024624



1094032



1025779



1020854

Description	Notes	Part no.
LCD computer monitor, 20"		1013805
LCD computer monitor, 24"	1900 x 1200 resolution	1043655
Computer, desktop for Alice systems		1043167
Computer, laptop for Alice systems		1024624
Speakers for Alice PSG		1094032
CPU/laptop microphone for Alice PSG		1025779
Broadway video card		1020854
Microphone kit, super-high gain		1097360


# Video/audio



1016347



1016348



1016349



1026576



1032204



1021731, 1021730



1077297

Description	Part no.
Video camera (fixed focus, analog)	1016348
Video web server axis model 2400+, 4 port	1016347
Analog video cable, 50 ft.	1021731
Analog video cable, 100 ft.	1021730
VideoEase CCTV balun	1026576
Modular CCTV power – through balun	1032204
Axis 215 PTZ network camera	1077297
IR projector	1016349
Microphone kit, super high-gain	1097360
Speakers for Alice PSG	1094032

# Printing/networking



9800



1012148

Description	Part no.
Printer, color laser, network	9800
16 port network switch	1012148

# Alice lab system accessories



1020266



1020584



1016343

Description	Notes	Part no.
Communication cable,Alice 6	Connects headbox to base station	1074689
CPAP therapy circuit interface kit,Alice 6		1082653
Headbox communication cable,Alice 6		1074689
Wall mounting kit,Alice 6		1082652
Shoulder strap for Alice 6 headbox		1078921
Wall mounting kit,Alice 5		1020266
Carrying pouch,Alice 5		1020584
Shoulder strap,Alice 5		1016343

# LoFlo capnography sensor



LoFlo connects to an Alice PSG system with a single USB connection and reduces the need to manage bulky equipment and calibration settings while patients wait to begin their sleep study.

Description	Part number
Alice, LoFlo C5 Sidestream module	1097955
Alice LoFlo adapter cable	1098078
Sampling line, 10 pk	3474-00
Alice, LoFlo starter kit, (LoFlo, sampling line 10 pk, Alice LoFlo USB cable, and cannulas)	1101087
Cannula, Pro-Flow, nasal, adult 10 pk	1097345
Cannula, Pro-Flow, nasal, pediatric 10 pk	1097348
Cannula, Pro-Flow, nasal, pediatric, small 10 pk	1097349
Cannula, Pro-Flow, nasal, pediatric 60 pk	1097939
Cannula, Pro-Flow, nasal, pediatric, small 60 pk	1097940

# Cables



H3150-50  
H3150-100



H3151-50  
H3151-100



1024583



8158-10



1000814



1016352



1016351



1007492



1021726, 1021727



8170-10



8192



8179




8155

Description	Notes	Part no.
Headbox communication cable, Alice 6		1074689
Connecting cable CAT-5 patch, 50'	Standard network patch cable	H3150-50
Connecting cable CAT-5 patch, 100'	Standard network patch cable	H3150-100
Connecting cable CAT-5 cross-over, 50'	Cable used for direct connections only, (i.e., PC to base station)	H3151-50
Connecting cable CAT-5 cross-over, 100'	Cable used for direct connections only, (i.e., PC to base station)	H3151-100
Patient cable, Alice 5	Connects base station to headbox	1024583
Sub-mini adapter for auxiliary input cable		8158-10
Encore shielded communications cable, 6'	Synchrony to Alice 5 base station	1000814
Connecting cable auxiliary input, 10'		1016352
A/B aux input signal splitter		1016351
SleepLink card with serial cable		1007492
3 ft. output analog cable	Connects Capnogard accessory devices to base station via analog connectors	1021726
6 ft. output analog cable	Connects Capnogard accessory devices to base station via analog connectors	1021727
Connecting cable to 8170 to 8130, 10'		8170-10
Parallel interface cable, 6'		8192
Reference electrode jumper cable		8179
Connecting cable 8130 auxiliary input to external device, 10'		8155
Connecting cable 8130 auxiliary input to external device, 25'		8155-25
Connecting cable 8130 auxiliary input to external device, 50'		8155-50
Connecting cable 8130 auxiliary input to external device, 75'		8155-75
Connecting cable 8130 auxiliary input to external device, 100'		8155-100



# Sensors – air flow



Description	Notes	Part no.
Thermistor airflow, adult, 1 channel, 249, Alice 6		1078754
Thermistor airflow, pediatric, 1 channel, 249, Alice 6		1078755
Thermistor airflow, pediatric small, 1 channel, 249, Alice 6		1078756
Thermistor, oral, adult, 249, Alice 6		1078757
Thermistor airflow, adult, 1 channel, 1.5mm, Alice 5		P1274
Thermistor airflow, pediatric, 1 channel, 1.5mm, Alice 5		P1336
Thermistor airflow, pediatric small, 1 channel, 1.5mm, Alice 5		P1337
Thermistor, oral, adult, 1.5mm, Alice 5		P1344
Pressure sensor, North America, Alice 6		1076590
Pressure sensor, International, Alice 6		1076592
PTAF Lite, Alice 5		P1304R
Nasal cannula, adult (Pro-Flow)	60 per pk	P1259
Nasal cannula, adult (Pro-Flow)	10 per pk	P1301
Nasal/oral cannula, adult (Pro-Flow Plus)	30 per pk	P1295
Nasal/oral cannula, adult (Pro-Flow Plus)	5 per pk	P1300
Nasal cannula, pediatric, small (Pro-Flow)	60 per pk	P1267
Nasal cannula, pediatric, small (Pro-Flow)	10 per pk	P1303
Nasal cannula, pediatric (Pro-Flow)	60 per pk	P1263
Nasal cannula, pediatric (Pro-Flow)	10 per pk	P1302

# Sensors – effort



1076333



P1736R



1073622



1075970



P1586, P1587



1094110



1094111

Description	Notes	Part no.
zRIP, DuraBelt kit, adult, Alice 6	Kit contains: 2 ea. P/N 1073622 zRIP DuraBelt adult effort belt and 2 ea. P/N P1133 wiresets	1078925
zRIP, DuraBelt kit, pediatric, Alice 6	Kit contains: 2 ea. P/N 1075970 zRIP DuraBelt pediatric effort belt and 2 ea. P/N P1133 wiresets	1078926
zRIP Inductance Effort, Alice 5/LE Adult	Kit contains: 2 ea. P/N 1073622 zRIP DuraBelt adult effort belts, P1736 zRIP Alice 5 driver module, 2 ea. P/N P1133 wiresets	1076333
zRIP Inductance Effort, Alice 5/LE Pediatric	Kit contains: 2 ea. P/N 1075970 zRIP DuraBelt pediatric effort belt, P1736 zRIP Alice 5 driver module, 2 ea. P/N P1133 wiresets	1076334
zRIP effort belt, adult		1073622
zRIP effort belt, pediatric		1075970
zRIP effort belt, large		1075973
zRIP module, Alice 5		P1736R
zRIP wireset	Qty: 1	P1133
ezRIP adult kit, Alice 5/LE, 84"	Kit contains: 1 abd ezRIP module, and 1 thx ezRIP module, 2 zRIP DuraBelts, adult	1099515
ezRIP module, abdomen, 1.5, 84", Alice 5/LE		1094110
ezRIP module, thorax, 1.5, 84", Alice 5/LE		1094111
CT2, adult, Alice 5	Piezo	P1586
CT2, pediatric, Alice 5	Piezo	P1587


# Sensors – oximetry



H5871



H5875



H5874



H5873



H5872



1019530, 1019531,  
1019532



1019529

Description	Notes	Part no.
<b>Alice 6</b>		
High performance/low power module		1075383
Masimo LNCS DCI adult reusable (clip)		1081984
Masimo LNCS DCIP, reusable (clip)		1081515
Masimo LNCS ADTX adhesive sensor, 20/box		1081517
Masimo LNCS PDTX adhesive sensor, 20/box		1081518
Masimo LNCS INF adhesive sensor, 20/box		1081519
Masimo LNCS YI, reusable		1081985
Masimo LNCS DBI, adult soft reusable		1081516
Masimo LNCS patient cable, LNC-4 ext		1082474
1597 Standard Multisite wrap, 100/box		1019531
1596 CleanShield Multisite wrap, 100/box		1019532
<b>Alice 5</b>		
Reusable oximeter sensor, finger clip		H5871
Disposable oximeter sensor, adult	Box of 20 ea.	H5875
Disposable oximeter sensor, pediatric	Box of 20 ea.	H5874
Disposable oximeter sensor, neonatal	Box of 20 ea.	H5873
Disposable oximeter sensor, preemie	Box of 20 ea.	H5872
Multisite oximeter probe, Masimo	Reusable probe, includes 6 ea. standard wraps	1019530
Multisite oximeter probe, Masimo	Standard wrap, box of 100 ea.	1019531
Multisite oximeter probe, Masimo	Clean shield wrap, box of 100 ea.	1019532
Patient cable, 8'	Use to interface all Masimo probes to Alice 5 headbox	1019529

## Sensors – body movement



P1694



1009963



1025751

Description	Notes	Part no.
Body position sensor, 249, Alice 6		1078758
Body position sensor, 1.5mm, Alice 5		P1694
PLM sensor, 249, Alice 6		1078752
Actimeter, Alice 5	Connects to the wrist, detects body movement	1009963
Actimeter tape	10 yds per roll	1025751

## Sensors – electrodes



1016360



1016361



1075397



1075402



8179



8051



8052

Description	Part no.
ECG, 72", DOM pack of 5 lead wires, domestic	1075232
ECG, 72", DOM pack of 5 lead wires, international	1075233
Adult repositionable ECG electrodes	1016360
Infant ECG electrodes (1 pkg. of 3 electrodes)	1016361
Reference electrode jumper cable, 10 ft.	8179
EEG electrode cream	8051
EEG electrode prep	8052
EEG/EMG/EOG, 48", gold cup	1075294
EEG/EMG/EOG, 72", gold cup	1075397
EEG/EMG/EOG, 96", gold cup	1075398
EEG/EMG/EOG, 120", gold cup	1075399
EMG, 72", snap, 10 per pack	1075400
EMG, 96", snap, 10 per pack	1075401
EMG, 120", snap, 10 per pack	1075402

# Sensors – snore



P1716, P1696, 1078753

Description	Part no.
Snore microphone, 249, Alice 6	1078753
Snore microphone, 1.5mm, Alice 5	P1716
Piezo snore sensor	P1696

## Miscellaneous accessories



1024582



1016309



1016310



1019279



1016344



1013223



8003-10



5963-00

Description	Notes	Part no.
Photo detector	For use with Alice 5 and Alice 6	1050795
Hard shell carrying case includes: sensor case		1024582
Sensor case, sensors not included		1016309
Mouse pad		1016310
Power supply, 60W, 6.25 VDC 3-6 AMPS, Alice 6		1019280
Uninterruptible power supply		1091779
Power supply 60W 6VDC, 10 AMP, Alice 5		1019279
Uninterrupted power supply (UPS)		1016344
Power cord, US domestic		1013223
Computer workstation, 36"		8003-10
Oxi-Pulse monitor	Includes power cord and manuals	1076305
Signal splitter Alice 5	For use with Oxi-Pulse	1016351
Single input AUX cable – 10 ft.	For use with Oxi-Pulse	1016352
Pro-Flow multi-purpose cannula, adult, 5 pk	For use with Oxi-Pulse	P1368
Pro-Flow multi-purpose cannula, ped, 5 pk	For use with Oxi-Pulse	P1369
Oxi-Pulse calibration kit	For use with Oxi-Pulse	1076306
Oxi-Pulse analog output cable	For use with Oxi-Pulse	1080848
Analog output module (AOM)		5963-00



# Software

Description	Part no.
Sleepware G3, English only	1082462
Sleepware G3, Multi-language	1082463
HL7 inbound only interface	1069193
HL7 outbound only interface	1069194
HL7 bi-directional interface	1067148

# Manuals

Description	Part no.
Sleepware G3 Getting Started Guide, English	1081398
Sleepware user guide, English	1016318
Sleepware user guide, English/International	1016319
Alice 6 hardware setup guide, English	1078160
Alice 6 hardware setup guide, English/International	1078161
Alice 5 hardware setup guide, English	1016329
Alice 5 hardware setup guide, English/International	1016330

# OmniLab Advanced + multi-mode titration system and accessories

With nine therapy modes, a quiet platform, customizable reports, user-friendly features and value-added benefits, you'll find the OmniLab Advanced + to be an ideal system for satisfying your lab titration needs. Along with OmniLab Direct and OmniLab Connect, the OmniLab Advanced + system provides a cohesive, harmonized solution to help you and your sleep technicians operate titration studies efficiently, easily and intuitively.



Description	Notes	Part no.
OmniLab Advanced + with Heated Tube humidifier		1111124
A Series System One humidifier		1111551
Water chamber (includes both 50 and 60 Series lids)		1063785
System One Heated Tube		SYS1HT15
OmniLab Connect		1111076
OmniLab Direct	Log in to <a href="http://My.Respironics.com">My.Respironics.com</a> or use PN to receive a CD	1112033

# Alice PDx


Alice PDx is a portable at-home diagnostic recording device for OSA screening, follow-up, and diagnostic assessment of cardio-respiratory sleep disorders. The system satisfies the AASM/CMS HST device guidelines and enables basic screening to comprehensive diagnostic evaluation. The Good Study Indicator conveys the amount of good quality data collected overnight to help you avoid or reduce the costs associated with retesting.



Alice PDx device kits (batteries, SpO<sub>2</sub> probe, and SD card reader, Sleepware also included)

Description	Part no.
Alice PDx device kit, (Domestic/US)	1043941
Alice PDx device kit, (International)	1043844
Alice PDx device kit, (Japan)	1047990
Alice PDx device kit, (China)	1047992

## Optional ExG expansion yoke



1040808, 1048415

1016339, 1020971

Description	Notes	Part no.
ExG yoke, Alice PDx (US and Japan)	Includes 1 ea. of P/N 1020971 and 1016339 and Sleepware analysis	1040808
ExG yoke, Alice PDx (International)	Includes 1 ea. of P/N 1020971 and 1016339 and Sleepware analysis	1048415
EEG/EMG/EOG leads 72 in., qty 10		1016339
EEG/EMG/EOG leads 48 in., qty 10		1020971

## Optional ECG expansion yoke



1040809, 1040810,  
1048916

1041199, 1041200

1016360

Description	Notes	Part no.
ECG yoke, Alice PDx (US)	Includes 1 ea. of P/N 1040809 and 5 ea. of P/N 1016360 and Sleepware analysis	1040809
ECG yoke, Alice PDx (International)	Includes 1 ea. of P/N 1041200 and 5 ea. of P/N 1016360 and Sleepware analysis	1040810
ECG yoke, Alice PDx (Japan)	Includes 1 ea. of P/N 1041200 and 5 ea. of P/N 1016360 and Sleepware analysis	1048916
ECG lead set 5-wire, 30 in. (76 cm), Domestic/US		1041199
ECG lead set 5-wire, 30 in. (76 cm), International		1041200
Adult ECG electrodes (pkg. of 5 electrodes)		1016360

# Sensors and replacement components



1075969



P1379



936



953



954A



P1388



P1328-60



P1343



1053279



1043937



1069860



1047300



1062049, 1062050,  
1062051



1065281

Description	Part no.
zRIP DuraBelt adult sensor kit, Alice PDx	1075969
Oral thermistor, Alice PDx	P1379
3 ft. SpO <sub>2</sub> extension cable (Nonin)	927
SpO <sub>2</sub> sensor (finger clip)	936
SpO <sub>2</sub> sensor (flex)	953
Adult flexwrap (disposable for use w/953 sensor, 25/pkg)	954A
zRIP DuraBelt effort sensor, adult single belt	1073622
zRIP wireset, Alice PDx	P1170
Thermistor, airflow, nasal/oral, adult, sensor, Alice PDx	P1388
Pro-Flow nasal cannula, adult, 16", 60 per case	P1328-60
Pro-Flow Plus nasal/oral cannula, adult, 16", 30 per case	P1343
CPAP titration kit, Alice PDx	P1391
Alice PDx LM sensor	Includes 1 ea. Alice PDx LM sensor, 1 Velcro strap PLM, small (20cm/18"), and 1 Velcro strap PLM, large (41cm/16")
Alice PDx battery door	1069622
Alice PDx holster	1053279
Alice PDx carry case	1053280
Secure digital card 1GB (2 SD cards per package)	1043937
Lanyard (replacement support lanyard)	1069860
SD card reader	1053194
Nonin finger soft sensor, large	1062049
Nonin finger soft sensor, medium	1062050
Nonin finger soft sensor, small	1062051
Surgical tape (box of 12 rolls of 10 yards ea.)	1065281

# Communication and therapy cables

Description	Part no.
Alice PDx sleep therapy cable	1040805
Alice PDx SleepLink cable	1040806
Alice PDx MobileLink cable	1040807
SleepLink card	1059963

## Manuals – professional



Description	Part no.
Alice PDx professional guide, (Domestic/US)	1059953
Alice PDx professional guide, (English International)	1050119
Alice PDx professional guide, (French)	1055341
Alice PDx professional guide, (Italian)	1055342
Alice PDx professional guide, (German)	1055343
Alice PDx professional guide, (Spanish)	1055344
Alice PDx professional guide, (Brazilian Portuguese)	1055345
Alice PDx professional guide, (Japanese)	1055346

## Manuals – user/patient



Description	Part no.
Alice PDx user/patient guide, (Domestic/US)	1059952
Alice PDx user/patient guide, (English International)	1050121
Alice PDx user/patient guide, (French)	1055347
Alice PDx user/patient guide, (Italian)	1055348
Alice PDx user/patient guide, (German)	1055349
Alice PDx user/patient guide, (Spanish)	1055350
Alice PDx user/patient guide, (Brazilian Portuguese)	1055351
Alice PDx user/patient guide, (Japanese)	1055352



# Educational materials

You can select from a variety of informative materials to help educate your staff and your patients.



Description	Part no.
Patient educational DVD	1055905
Alice PDx interactive tutorial for clinicians	1055918
Alice PDx case overlay (4/pkg)	1053948
Alice PDx patient instruction reference card	1065860
What you need to know about RIP poster	1081437
zRIP FAQ learning guide	1080365

# Stardust II

This portable at-home diagnostic system is powerful enough to generate comprehensive seven-channel reports, yet small enough to be worn comfortably while sleeping – wherever, whenever. The Stardust II also enables you to meet level III and IV testing requirements.



Stardust II

Description	Part no.
Stardust II recorder package, US	1025720
Stardust II recorder package, International	1011176
Stardust II with zRIP DuraBelt, US	1098874

## Cables



1006037



5910



612214K



1007492



1011955

Description	Part no.
Optical remote cable, 15"	1006037
Tranquility Auto interface	5910
Serial connection RS 232 download cable	612214K
SleepLink card with serial cable	1007492
REMstar auto serial link cable	1011955

## Sensors



1099558



1094097



P1805

Effort			
Description	Notes	Part no.	
ezRIP adult, Stardust II	Kit contains: 1 thx ezRIP module, 1 zRIP DuraBelt, adult	1099558	
ezRIP module, thorax, harwin, 26"		1094097	
CT2 effort sensor, adult, Stardust		P1805	

# Sensors



1007199



612054



1101838

Oximetry		
Description	Notes	Part no.
Oximeter finger clip sensor	Replaces Y sensor	1007199
Flexible finger oximetry Y probe kit	Includes adult sensor tape	612054
Butterfly wrap, adult	10 per pk	1101838



P1396



P1340



P1339



P1328-60



P1343

Airflow		
Description	Notes	Part no.
Stardust PTAF sensor kit	Includes: 1 sensor and 5 cannulas	P1396
Stardust PTAF sensor		P1340
Thermistor sensor, adult, Stardust		P1339
Pro-Flow adult nasal cannula, 16"	60 per pk	P1328-60
Pro-Flow Plus adult nasal/oral cannula, 16"	30 per pk	P1343

## Miscellaneous



1008312



1003510



1063504

Description	Part no.
Pouch	1008312
Stardust II carrying case	1064344
Respironics screwdriver	1003510
Stardust II patient instruction DVD	1063504

## Manuals

Description	Part no.
Stardust II recorder manual	1026572
Host software manual v2.0, Domestic	1026567
Host software manual v2.0, International English	1026566
Host software manual v2.0, French	1018899
Host software manual v2.0, German	1026564
Host software manual v2.0, Spanish	1026565

# Actigraphy systems

Our Actiwatch Spectrum and Actiwatch 2 ambulatory products are designed to help make actigraphy technology easier to use, the devices more comfortable to wear, and the analysis functions less complex and time consuming.



Description	Notes	Part no.
Actiwatch 2 startup system	Includes device, dock system, software, watchband kit, clinician, technician and wearer guides	1044809
Actiwatch Spectrum startup system	Includes device, dock system, software, watchband kit, clinician, technician and wearer guides	1046964
Actiwatch 2	Includes device and wearer guide	1048090
Actiwatch 2 dock communication system	Includes dock, USB cable, power adapter and supply, tech guide, does not include the Actiwatch 2	1048092
Actiwatch Spectrum	Includes device and wearer guide	1048091
Actiwatch Spectrum dock communication system	Includes dock, USB cable, power adapter and supply, tech guide, does not include the Actiwatch Spectrum	1048093
Actiwatch Spectrum PRO	Includes device and wearer guide	1099351
Actiwatch Spectrum Plus	Includes device and wearer guide	1101894
Actiwatch Spectrum PRO startup system	Includes device, USB cable, software, watchband kit, clinician, technician and wearer guides	1109515
Actiwatch Spectrum Plus startup system	Includes device, USB cable, software, watchband kit, clinician, technician and wearer guides	1109516
Actiware software		1104775
Actiware additional site license		1114828

# Actigraphy accessories

Description	Notes	Part no.
Actiwatch 2 watchband kit	Includes tool, 3 watchbands, 6 pins and tool kit	1048135
Actiwatch Spectrum watchband kit	Includes tool, 3 watchbands, 6 pins and tool kit	1048135
Watchband tool		1044811
Actiwatch 2 watchband pack	Includes 3 watchbands and 6 pins	1044810
Actiwatch Spectrum watchband pack	Includes 3 watchbands and 6 pins	1048136
Actiwatch 2 hospital band pack	Includes 25 hospital bands and 12 spring bars	1044816
Actiwatch power adapter		1044812
Actiwatch Spectrum and Actiwatch 2 USB cable		1044813
Actiwatch 2 power supply		1044814
Actiwatch 2 docking station		1046977
Actiwatch Spectrum Plus/PRO USB cable	Type A to micro B	1104765





# Pro-tech sensors

As a leader in the development, manufacturing, and supply of sleep diagnostic sensors, we offer a full line of high-performance, quality Pro-Tech sensors to help you maintain low cost of ownership while obtaining the data you need. Our sensor product line includes durable and reliable sensors that detect respiratory airflow, respiratory effort, snoring, periodic limb movements, and sleep position.

Not all products are available in all countries. Check with Philips Respironics customer service or your local distributor.

**Pro-Flow nasal cannulas**

Description	Part number
Pro-Flow, nasal, adult, 7', (60 per case), male luer lok	P1259
Pro-Flow, nasal, adult, 7', (10 pack), male luer lok	P1301
Pro-Flow, nasal, adult, 8.25", (60 per case), male luer lok	P1306
Pro-Flow, nasal, adult, 1.1M, (20 per case), male luer lok	P1317-20
Pro-Flow, nasal, adult, 16", (60 per case), male luer lok	P1328-60
Pro-Flow, nasal, adult, 16", (5 pack), male luer lok	P1390
Pro-Flow, nasal, adult, 7', (10 pack), no filter, male luer lok	1091077
Pro-Flow, nasal, adult, 7', (60 per case), no filter, male luer lok	1091078
Pro-Flow, nasal, pediatric, 7', (60 per case), male luer lok	P1263
Pro-Flow, nasal, pediatric, 7', (10 pack), male luer lok	P1302
Pro-Flow, nasal, pediatric, 1.1M, (20 per case), male luer lok	P1321-20
Pro-Flow, nasal, pediatric, 16", (60 per case), male luer lok	P1381
Pro-Flow, nasal, pediatric small, 7', (60 per case), male luer lok	P1267
Pro-Flow, nasal, pediatric small, 7', (10 pack), male luer lok	P1303
Pro-Flow, nasal, pediatric small, 1.1M, (20 per case), male luer lok	P1323-20



**Pro-Flow and Pro-Flow Plus cannulas**

- Designed to enhance nasal and oral signal transfer
- Soft and narrow nasal prongs help reduce nasal resistance and increase patient comfort
- Incorporate narrow kink-resistant tubing to help facilitate excellent response to high-frequency signals such as snoring
- Hydrophobic, anti-microbial inline filter
- Pediatric and adult sizes

**Pro-Flow Plus nasal/oral cannulas**

Description	Part number
Pro-Flow Plus, nasal/oral, adult, 7', (30 per case), male luer lok	P1295
Pro-Flow Plus, nasal/oral, adult, 7', (5 pack), male luer lok	P1300
Pro-Flow Plus, nasal/oral, adult, 8.25", (30 per case), male luer lok	P1308
Pro-Flow Plus, nasal/oral, adult, 8.25", (30 per case), no filter, male luer lok	P1334
Pro-Flow Plus, nasal/oral, adult, 1.1M, (10 per case), male luer lok	P1319-10
Pro-Flow Plus, nasal/oral, adult, 16", (30 per case), male luer lok	P1343
Pro-Flow Plus, nasal/oral, pediatric, 7', (30 per case), male luer lok	P1354
Pro-Flow Plus, nasal/oral, pediatric, 16", (30 per case), male luer lok	P1383



**Pro-Flow Plus nasal/oral cannulas**

- Addresses monitoring of mouth-breathing patients
- Oral prong detects pressure changes in the mouth
- Pediatric and adult sizes

## Pressure airflow (continued)

### Pro-Flow MultiPurpose nasal cannulas

Description	Part number
Pro-Flow MultiPurpose, nasal, adult, 7', (30 per case), male luer lok	P1350
Pro-Flow MultiPurpose, nasal, adult, 7', (5 pack), male luer lok	P1368
Pro-Flow MultiPurpose, nasal, pediatric, 7', (30 per case), male luer lok	P1346
Pro-Flow MultiPurpose, nasal, pediatric, 7', (5 pack), male luer lok	P1369
Pro-Flow MultiPurpose, nasal, pediatric small, 7', (30 per case), male luer lok	P1348
Pro-Flow MultiPurpose, nasal, pediatric small, 7', (5 pack), male luer lok	P1370



### Pro-Flow Plus MultiPurpose cannulas

Description	Part number
Pro-Flow Plus MultiPurpose, nasal/oral, adult, 7', (30 per case), male luer lok	P1373
Pro-Flow Plus MultiPurpose, nasal/oral, pediatric, 7', (30 per case), male luer lok	P1366



### Cannula accessories

Description	Part number
Cannula adapter kit, female luer thread style on 1/8 inch tubing	P1324
Double-port CPAP titration kit Kit contains: 12 cannula titration filters, 2 male luer loks, 2 circuit connection ports, small bore connecting tubing	P1327
CPAP titration kit, Alice PDx, connects CPAP circuit to PDx via male luer lok	P1391
Cannula or titration filter, (50 pack), male and female luer lok	P1107
Cannula or titration filter, (10 pack), male and female luer lok	P1109
RUSleeping RTS kit, male luer lok	1041483
Pro-Wrap, 100 Pro-Wrap per roll, no active sensor	P1162
Pro-Flow sample pack, includes P1259 (qty: 2) and P1295 (qty: 2)	P1258

### Pro-Flow MultiPurpose and Pro-Flow Plus MultiPurpose cannulas

- Designed to concurrently monitor airflow and end-tidal CO<sub>2</sub>
- Provides ability to monitor three parameters (airflow, snore, and EtCO<sub>2</sub>) with one patient interface
- Can be used with most pressure transducers and continuous sampling capnographs
- Hydrophobic, anti-microbial inline filter
- Pediatric and adult sizes

# Thermal airflow

## Reusable thermocouple (cannula style)

Description	Part number
Thermocouple 6" single element, adult, 1.5mm connectors	P1221
Thermocouple nasal/oral airflow sensor, 1 channel, adult, 1.5mm connectors	P1222
Thermocouple nasal/oral airflow sensor, 2 channel, adult, 1.5mm connectors	P1224
Thermocouple nasal/oral airflow sensor, 3 channel, adult, 1.5mm connectors	P1226
Thermocouple oral airflow sensor, adult, 1.5mm connector	P1273
Thermocouple airflow sensor, 1 channel, adult, binder connector	P1757
Thermocouple airflow sensor, 1 channel, adult, binder connector	P1758
Thermocouple nasal/oral airflow sensor, 1 channel, pediatric, 1.5mm connectors	P1223
Thermocouple nasal/oral airflow sensor, 2 channel, pediatric, 1.5mm connectors	P1225
Thermocouple airflow sensor, 1 channel, pediatric, binder connector	P1759
Thermocouple airflow sensor, 1 channel, pediatric, binder connector	P1760
Thermocouple nasal airflow sensor, 1 channel, pediatric small, nasal bridgemount, 1.5mm connectors	P1214
Thermocouple airflow sensor, 1 channel, pediatric small, binder connector	P1761
Thermocouple airflow sensor, 1 channel, pediatric small, binder connector	P1762



### Thermocouples and thermistors

- Compatible with most major PSG systems
- Ability to monitor oral/nasal separately and/or together with adult and pediatric sizes and disposable thermocouples
- Insulated tips
- Six-month warranty

## Disposable thermocouple

Description	Part number
Reusable wireset	P1238
Disposable thermocouple nasal/oral airflow sensor, adult, (10 pack), 1.5mm connectors	P1242
BreathSensor cable, adult, 1.5mm connectors	P1332
BreathSensor cable, adult, 249 connector (2-pin keyhole) (orange)	P1154
BreathSensor cable, SD20, adult, 220 connector (3-pin keyhole)	P1112

## Thermal airflow (continued)

### Reusable thermistor

Description	Part number
<b>Alice 3</b>	
Thermistor; 1 channel, nasal/oral, adult, 2.5 connector	P1246
Thermistor 1 channel, nasal/oral, pediatric, 2.5mm connector	P1298
Thermistor 1 channel, pediatric small, nasal bridgemount, 2.5mm connector	P1296
<b>Alice 4</b>	
Thermistor; 1 channel, nasal/oral, adult, LEMO triaxial connector	P1286
Thermistor 1 channel, nasal/oral, pediatric, LEMO triaxial connector	P1299
Thermistor 1 channel, pediatric small, nasal bridgemount, LEMO triaxial connector	P1297
<b>Alice 5</b>	
Thermistor; 1 channel, nasal/oral, adult, 1.5mm connectors	P1274
Thermistor 1 channel, nasal/oral, pediatric, 1.5mm connectors	P1336
Thermistor 1 channel, pediatric small, nasal bridgemount, 1.5mm connectors	P1337
Thermistor; oral, adult, 1.5mm connectors	P1344
<b>Alice 6</b>	
Thermistor; 1 channel, nasal/oral, adult, 249 connector (2-pin keyhole)	1078754
Thermistor; 1 channel, nasal/oral, pediatric, 249 connector (2-pin keyhole)	1078755
Thermistor 1 channel, pediatric small, nasal bridgemount, 249 connector (2-pin keyhole)	1078756
Thermistor; oral, adult, 249 connector (2-pin keyhole)	1078757
<b>Alice PDx</b>	
Thermistor airflow, nasal/oral, adult, Alice PDx, 2.5mm connector	P1388
Thermistor; oral, adult, Alice PDx, 2.5mm connector	P1379
<b>Stardust II</b>	
Thermistor airflow, nasal/oral, adult, Stardust II, harwin connector	P1339
<b>Other</b>	
Thermistor; 1 channel, nasal/oral, adult, Sandman SD20, 220 connector (black)	P1364
Thermistor; 1 channel, nasal/oral, pediatric, Sandman SD20, 220 connector (black)	P1378
Thermistor; 1 channel, nasal/oral, adult, Embletta, 249 connector (2-pin keyhole) (gray)	P1274-249
Thermistor; 1 channel, nasal/oral, pediatric, Embletta, 249 connector (2-pin keyhole) (gray)	P1298-249



## Pressure transducers

### PTAF2

Description	Part number
PTAF2 (kit contains: PTAF2 module, a 9V battery, and one P1258 cannula sample pack)	P1287

### PTAF2 cables

Description	Part number
High level cable, 7', 1/8 male connector	P1078
High level cable, 7', 1/8 female connector	P1079
High level cable, 7', 3/32 male connector	P1080
High level cable, 7', 1/4 male connector	P1081
High level cable, 7', 6-4 modular connector	P1082
High level cable, 7', 3.5mm connector	P1138
Low level cable, 7', 1/8 male connector	P1076
Low level cable, 7', 3/32 male connector	P1077
Low level cable, 7', 1.5mm connector	P1075
Extension with coupler, 25'	P1087-025
Extension with coupler, 50'	P1087-050
Extension with coupler, 100'	P1087-100



### PTAF2 pressure-based transducer

- Designed to generate accurate respiratory signals
- Includes high-level and low-level outputs, unfiltered or filtered airflow output, snoring output, and remote-control options
- Used in AC or DC configurations and during CPAP titrations
- One-year warranty



## Pressure transducers (continued)

### Pressure sensor

Description	Part number
Pressure sensor; domestic, Alice 6	1076590
Pressure sensor; international, Alice 6	1076592



### PTAF Lite

Description	Part number
PTAF Lite startup kit Kit contains: P1258: cannula sample pack; P1289: PTAF Lite	P1289
PTAF Lite for Alice 5	P1304R



### PTAF2 and PTAF Lite accessories

Description	Part number
Cannula titration filter (50 pack)	P1107
Cannula titration filter (10 pack)	P1109
Cannula adapter kit	P1324
CPAP titration kit	P1327
CPAP titration kit Alice PDx	P1391
Cannula sample pack	P1258

### Stardust PTAF

Description	Part number
Stardust PTAF sensor kit, (includes 5 cannula starter pack), harwin connector	P1396
Stardust PTAF sensor (sensor only), harwin connector	P1340

The PTAF Lite is a small, lightweight piezo-style pressure transducer used in AC device configurations.

Please have available the make and model number of polygraph you are using when you contact our customer service department to receive the correct configuration for your unit.

## Respiratory inductance plethysmography effort sensors

### zRIP DuraBelt startup kits

Description	Part number
<b>Alice 3</b>	
<b>zRIP DuraBelt kit, adult, Alice 3</b> Kit contains: 1 each P1734 zRIP Alice 3 driver module, 2 each P1133 zRIP wiresets, 2 each 1073622 zRIP DuraBelts, adult	1076331
<b>zRIP DuraBelt kit, pediatric, Alice 3</b> Kit contains: 1 each P1734 zRIP Alice 3 driver module, 2 each P1133 zRIP wiresets, 2 each 1075970 zRIP DuraBelts, pediatric	1076332
<b>Alice 5</b>	
<b>zRIP DuraBelt kit, adult, Alice 5</b> Kit contains: 1 each P1736R zRIP Alice 5, 2 each P1133 zRIP wiresets, 2 each 1073622 zRIP DuraBelts, adult	1076333
<b>zRIP DuraBelt kit, pediatric, Alice 5</b> Kit contains: 1 each P1736R zRIP driver module, 2 each P1133 zRIP wiresets, 2 each 1075970 zRIP DuraBelts, pediatric	1076334
<b>Alice 6</b>	
<b>zRIP DuraBelt kit, adult, Alice 6</b> Kit contains: 2 each 1073622 zRIP DuraBelts, adult, 2 each P1133 zRIP wiresets	1078925
<b>zRIP DuraBelt kit, pediatric, Alice 6</b> Kit contains: 2 each 1075970 zRIP DuraBelts, pediatric, 2 each P1133 zRIP wiresets	1078926
<b>Alice PDx</b>	
<b>zRIP DuraBelt kit, adult, Alice PDx</b> Kit contains: 2 each P1170 wiresets, 2 each 1073622 zRIP DuraBelt, adult	1075969
<b>Other</b>	
<b>zRIP DuraBelt kit, adult</b> Kit contains: 1 each P1732 zRIP driver module, 2 each P1133 zRIP wiresets, 2 each 1073622 zRIP DuraBelt, adult	1076329
<b>zRIP DuraBelt kit, pediatric</b> Kit contains: 1 each P1732 zRIP driver module, 2 each P1133 zRIP wiresets, 2 each 1075970 zRIP DuraBelts, pediatric	1076330
<b>zRIP DuraBelt Sum RT kit, adult</b> Kit contains: 1 each P1753 zRIP Sum RT driver module, 2 each P1133 zRIP wiresets, 2 each 1073622 zRIP DuraBelts, adult	1076336
<b>zRIP DuraBelt Sum RT kit, pediatric</b> Kit contains: 1 each P1753 zRIP driver module, 2 each P1133 zRIP wiresets, 2 each 1075970 zRIP DuraBelts, pediatric	1076337
<b>zRIP DuraBelt kit, adult, DC</b> Kit contains: 1 each P1738 zRIP driver module, 2 each P1133 zRIP wiresets, 2 each 1073622 zRIP DuraBelts, adult	1076335



### zRIP respiratory effort sensors and zRIP DuraBelts

- Provide linear breathing movement tracings with no false paradoxing
- Incorporate enhanced woven wire and buckle design
- Offer complete patient coverage
- Facilitate easy initial placement, and easy middle-of-the-night adjustments
- Fully washable
- Available in non-SUM and SUM RT (breath-to-breath summation) modules
- Replaceable batteries
- One-year warranty

## Respiratory inductance plethysmography effort sensors (continued)

### zRIP DuraBelt startup kits (continued)

Description	Part number
<b>Other (continued)</b>	
<b>zRIP DuraBelt kit, pediatric, DC</b> Kit contains: 1 each P1738 zRIP driver module, 2 each P1133 zRIP wiresets, 2 each 1075970 zRIP DuraBelts, pediatric	1075972
<b>zRIP DuraBelt kit, adult, QDC</b> Kit contains: 1 each P1176 abdomen cable, 1 each P1177 ribcage cable and 2 each 1073622 zRIP DuraBelts, adult	1075968
<b>zRIP DuraBelt kit, pediatric, QDC</b> Kit contains: 1 each P1176 abdomen cable, 1 each P1177 ribcage cable and 2 each 1075970 zRIP DuraBelts, pediatric	1075967



### zRIP modules

Description	Part number
zRIP module, 1.5mm connectors	P1732
zRIP SUM RT module, 1.5mm connectors	P1753
zRIP module, Alice 5, 1.5mm connectors	P1736R
zRIP module, Alice 3, LEMO, 2-pin connector	P1734
zRIP module, DC, 3.5mm connector	P1738



### zRIP DuraBelt effort sensors

Description	Part number
zRIP DuraBelt adult, size: 24" - 74", 1.5mm panel mount black	1073622
zRIP DuraBelt adult, large, size: 45" - 123", 1.5mm panel mount black	1075973
zRIP DuraBelt, pediatric, size: 16" - 42", 1.5mm panel mount black	1075970

### zRIP wiresets and cables

Description	Part number
zRIP wireset, 249 connector (2-pin keyhole)	P1133
zRIP wireset, Medicare, 249 connector (2-pin keyhole)	P1137
zRIP wireset, Alice PDx, 249 connector (2-pin keyhole)	P1170
zRIP snap adapter	P1151
zRIP cable abdomen, proprietary QDC connector	P1176
zRIP cable ribcage, proprietary QDC connector	P1177

## Respiratory inductance plethysmography effort sensors (continued)

### ezRIP DuraBelt kits

Description	Part number
<b>ezRIP adult kit, 1.5mm, 84"</b> 1 Abd, 1 Thx, 2 zRIP DuraBelts, adult	1099514
<b>ezRIP adult kit, Alice 5/LE, 84"</b> 1 Abd, 1 Thx, 2 zRIP DuraBelts, adult	1099515
<b>ezRIP adult kit, 249, 84"</b> 1 Abd, 1 Thx, 2 zRIP DuraBelts, adult	1099516
<b>ezRIP adult kit, 249, 16"/20"</b> 1 Abd, 1 Thx, 2 zRIP DuraBelts, adult	1099517
<b>ezRIP adult, Stardust II kit</b> 1 Thx, 1 zRIP DuraBelt, adult	1099558

### ezRIP combination module/wireset (individual component)

Description	Part number
ezRIP module, abdomen, 1.5mm, 84"	1094094
ezRIP module, thorax, 1.5mm, 84"	1094096
ezRIP module, abdomen, 1.5, 84", Alice 5/LE	1094110
ezRIP module, thorax, 1.5, 84", Alice 5/LE	1094111
ezRIP module, abdomen, 249, 84"	1094093
ezRIP module, thorax, 249, 84"	1094095
ezRIP module, abdomen, 249, 20"	1094920
ezRIP module, thorax, 249, 16"	1094921
ezRIP module, thorax, harwin, 26"	1094097



### ezRIP combination module/wireset sensor

- Low initial RIP investment
- Minimizes bedside module congestion
- Easy configurability – the driver module clips directly onto the zRIP DuraBelt and legacy zRIP belt
- Non-replaceable batteries
- One-year warranty

### Understanding RIP signals vs. piezo signals

RIP	Piezo
Directly reflects breathing activity of a patient	Indirectly reflects breathing activity of a patient
Sensor wire wraps around the patient	Sensor (piezo) most commonly located near the belt buckle
Belt only needs to be snug on the patient for an accurate signal	Belt needs to be pulled taut for an accurate signal
Gain settings do not change during the night when the signal decreases	Gain settings usually need to be altered during the night to capture the signal
The signal decreasing reflects true patient breathing effort	The signal decreasing, or changing, can reflect “trapping” or artificial artifact
Thoracic belt placement is designed to be more precise than piezo	Belt placement can be anywhere in the thoracic area and pulled tight

## Piezo respiratory effort sensors

### Piezo respiratory effort kits

Description	Part number
<b>Piezo 1 channel adult kit</b> Kit contains: 1 each P1420 adult Velcro® strap kit 1.5", 1 each P1460 piezo sensor 1.5", 1.5mm connectors	P1461R
<b>Piezo 2 channel plus adult kit</b> Kit contains: 2 each P1420 adult Velcro strap kit 1.5", 3 each P1460 piezo sensor 1.5", 1.5mm connectors	P1463R
<b>Piezo 1 channel pediatric kit</b> Kit contains: 1 each P1416 pediatric Velcro strap kit 1.5", 1 each P1460 piezo sensor 1.5", 1.5mm connectors	P1483
<b>Piezo 1 channel pediatric kit</b> Kit contains: 1 each P1448 pediatric Velcro strap kit 1", 1 each P1453 pediatric piezo sensor 1", 1.5mm connectors	P1462R
<b>Piezo 2 channel plus pediatric kit</b> Kit contains: 2 each P1448 pediatric Velcro strap kit 1", 3 each P1453 pediatric piezo sensor 1", 1.5mm connectors	P1464R
<b>Adult Velcro strap kit, 1.5"</b> Kit contains: 1 each adult Velcro strap 1.5", size 1, size 2, and size 3	P1420
<b>Pediatric Velcro strap kit, 1.5"</b> Kit contains: 1 each pediatric Velcro strap 1.5", size 1, size 2, and size 3	P1416
<b>Pediatric Velcro strap kit 1"</b> Kit contains: 1 each pediatric Velcro strap 1", size 1, size 2, size 3, and size 4	P1448

### Piezo respiratory effort sensors

Description	Part number
Piezo sensor, 1.5", adult, 1.5mm connectors	P1460
Piezo sensor, 1" pediatric, 1.5mm connectors	P1453
Piezo sensor, 1.5" adult, Alice 4, LEMO triaxial connector	P1571
Piezo sensor, 1.5" adult, Alice 3, 2-pin LEMO connector	P1547
CT2 piezo sensor, adult, 1.5mm connectors	P1582
CT2 piezo sensor, pediatric, 1.5mm connectors	P1585
Stardust CT2 effort sensor, adult, harwin connector	P1805
CT2 effort sensor, adult, Alice 5, 1.5mm connectors	P1586
CT2, effort sensor, pediatric, Alice 5, 1.5mm connectors	P1587
CT2, effort sensor, adult, Alice 4, LEMO triaxial connector	P1804
CT2, effort sensor, pediatric, Alice 4, LEMO triaxial connector	P1814
CT2, effort sensor, adult, Alice 3, 2-pin LEMO connector	P1803
CT2, effort sensor, pediatric, Alice 3, 2-pin LEMO connector	P1813
CT2, effort sensor, adult, Embletta, 249 connector (2-pin keyhole)	P1800
CT2, effort sensor, chest, adult, dark gray, 16", 249 connector (2-pin keyhole)	P1823
CT2, effort sensor, abdominal, adult, light gray, 19", 249 connector (2-pin keyhole)	P1824
CT2, effort sensor, chest, pediatric, dark gray, 16", 249 connector (2-pin keyhole)	P1825
CT2, effort sensor, abdominal, pediatric, light gray, 19", 249 connector (2-pin keyhole)	P1826



### Piezo respiratory effort sensors

- Sensor straps are easy to size on most patients
- Low profile of the sensor straps enhances patient comfort
- Compatible with most PSG systems
- Plugs directly into standard AC channels
- Latex free
- One-year warranty



### Piezo respiratory effort accessories

Description	Part number
CT2, 18" extension strap	P1592
CT1 belt strap, adult, size 1, 1.5"w, 33.5"l	P1421
CT1 belt strap, adult, size 2, 1.5"w, 49"l	P1422
CT1 belt strap, adult, size 3, 1.5"w, 77"l	P1423

## Body position sensors

### Body position sensors

Description	Part number
SPI (sleep position indicator) resistive, Cadwell, 8-8 modular plug	P1652
Body position sensor, DC with 3.5mm connector	P1664
Body position sensor, DC with 1.5mm connector	P1664-WS2074
SPI resistive Flaga, 1.5mm connector	P1676
Body position sensor, 220 connector (3 pin keyhole), Bio-Logic	P1692
Body position sensor, Alice 6, 249 connector (2-pin keyhole)	1078758
Body position sensor, Alice 5, 1.5mm connectors	P1694
Body position sensor, Alice 4, LEMO triaxial connector	P1707
Body position sensor, Alice 3, LEMO connector	P1706
Body position sensor, Sandman SD20, 249 connector (2-pin keyhole) (yellow)	P1719
Body position sensor, Sandman Pocket, 249 connector (2-pin keyhole) (gray)	P1744
Body position sensor, DC with 3.5mm female connector	P1746



### Body position sensors

- Designed with patient comfort in mind
- Ability to identify the primary sleep positions (prone, supine, left, right, and upright)
- Patient's body position can be continuously monitored regardless of the amount of movement
- Attaches to the thoracic respiratory effort belt easily and comfortably
- Compatible with most PSG systems
- One-year warranty



## Snore sensors and microphones

### Snore sensors and microphones

Description	Part number
Piezo snore sensor, 1.5mm connector	P1696
Snore sensor, Embletta, 249 connector (2-pin keyhole)	P1718
Snore microphone, 1.5mm connector	P1716
Snore microphone, Alice 6, 249 connector (2-pin keyhole)	1078753
Snore microphone, Alice 3, 3/32 phone connector	8017
Microphone snore, Sandman Pocket, 249 connector (2-pin keyhole)	P1742
Microphone snore, binder, 50cm wire length	P1900



### Snore sensors and microphones

- The snore sensor, with its thin piezo crystal, is designed to provide an accurate, reproducible signal
- The snore microphone is designed to provide a thin baseline and clear, responsive signals
- Compatible with most PSG systems
- Six-month warranty for snore sensors
- Ninety-day warranty for snore microphone

## Periodic limb movement

### Periodic limb movement (PLM)

Description	Part number
PLM sensor, Alice 6, 249 connector (2-pin keyhole)	1078752
Piezo PLM, 1.5mm connectors	P1697
<b>Piezo PLM kit</b>	P1698R
Kit contains: 1 each P1697 PLM sensor, 1 each P1689 Velcro strap PLM large, 1 each P1690 Velcro strap PLM small	
<b>Piezo PLM 2 per kit</b>	P1699R
Kit contains: 2 each P1697 PLM sensor, 2 each P1689 Velcro strap PLM large, 2 each P1690 Velcro strap PLM small, 1 each, P1097 safety male jumper wireset	

### Periodic limb movement (PLM) accessories

Description	Part number
Jumper, wireset, 1.5mm	P1097
Velcro strap PLM, large	P1689
Velcro strap PLM, small	P1690



### Periodic limb movement sensors

- Detect a patient's leg or arm movement
- Designed to provide clean, reliable signals while providing exceptional patient comfort
- Sensitive to subtle movements
- Do not require skin preparation prior to application
- Using kits, multiple sensors can be configured to combine PLM signals onto a single channel
- Six-month warranty

# Sensor tester

## Sensor tester

Description	Part number
ST-3, tests most PSG, includes P1315, continuity test wireset	P1711R
ST-3, Alice 5, includes P1315, continuity test wireset	P1722

## Sensor tester accessories

Description	Part number
ST-3, continuity test wireset	P1315
Low level cable, 1.5mm connectors	P1075
High level cable, 1/8 male connector	P1078
High level cable, 3.5mm connector	P1138



### ST-3 sensor tester

- Provides ability to test sensors prior to sleep studies
- Checks and verifies EEG electrodes and wires for continuity
- Square wave generator verifies polygraph channels and helps verify that the jackbox inputs and DC hi-level boxes are functional
- One-year warranty

## Education materials

Beyond just offering products and technology, we help you help staff and patients by providing product literature, instructional videos, reimbursement information, and patient educational materials.

Category	Description	Format	Part number
<b>For staff</b>			
Product support	Alice PDx clinician tutorial	CD-ROM	1055918
	What you need to know about RIP	Poster	1081437
	Actigraphy clinical scenarios	DVD	1043293
	RIP and zRIP learning guide and FAQ's	Booklet	1079901
Reimbursement information	Helpful hints for filing – Home sleep testing sleep studies	Brochure	1075752
	Helpful hints for filing – Polysomnography and home sleep testing for diagnosing OSA	Brochure	1093107
	Helpful hints for filing – Oral appliances for the treatment of OSA	Brochure	1091987
Continuing education	Intro to sleep diagnostic equipment monograph	Booklet	1052673
Insomnia program	Cognitive behavioral therapy for insomnia (CBTI)	Binder	Ask your sales representative
<b>For patients</b>			
Patient education	Having a sleep study	Pamphlet	1049375
	Having a home sleep study	Pamphlet	1064487
	I have sleep apnea. Now what?	Brochure	1038370
Product setup	Stardust II patient instructions, video	DVD	1063504
Product setup	Alice PDx patient instructions, video	DVD	1055905



## Appendix: AASM scoring manual guidelines and Pro-Tech sensor application

AASM scoring manual guideline – adults	Pro-Tech sensor
<b>Thermal airflow</b> Adults VIII.1.A: “The sensor to detect absence of airflow...is an oronasal thermal sensor.”	Oronasal thermocouple P1222
<b>Flow nasal pressure</b> Adults VIII.1.B: “The sensor for detection of... hypopnea is a nasal air pressure transducer...”	PTAF2 startup kit P1287 Nasal cannulas: adult (60 per case) P1259
<b>Respiratory effort</b> Adults VIII.1.C: “The sensor for detection of respiratory effort is.... inductance plethysmography.”	Respiratory Inductance Plethysmography: zRIP DuraBelt kit, adult 1076329 or ezRIP adult kit, 1099514

AASM scoring manual guideline – children	Pro-Tech sensor
<b>Thermal airflow</b> Children VIII.1.A: “The sensor to detect absence of airflow...is an oronasal thermal sensor.”	Oronasal thermocouple P1223
<b>Flow nasal pressure</b> Children VIII.1.B: “The sensor for detection of... hypopnea is a nasal air pressure transducer...”	PTAF2 startup kit P1287 Nasal cannulas: pediatric (60 per case) P1263
<b>Respiratory effort</b> Children VIII.1.C: “The sensor for detection of respiratory effort is.... inductance plethysmography.”	Respiratory Inductance Plethysmography: zRIP DuraBelt kit, pediatric 1076330 or ezRIP module, abdomen, 1094094, ezRIP module, thorax, 1094096 and zRIP DuraBelt, pediatric, 1075970
<b>Hypoventilation</b> Children VIII.1.E: “Acceptable methods for assessing alveolar hypoventilation.... end-tidal PCO <sub>2</sub> monitoring.” Enable EtCO <sub>2</sub> and pressure monitoring simultaneously from both nostrils using one cannula.	Pro-Flow MultiPurpose Cannula: pediatric (30 per case) P1346





# Value-added programs

Philips Respironics offers a broad array of programs aimed at helping to maximize your growth and operational efficiencies. From compliance and resupply to training and warranties, we provide a spectrum of unique programs focused on facilitating success for you and your patients.

Not all products are available in all countries. Check with Philips Respironics customer service or your local distributor.



## Device programs

### Flex Promise Program

Our Flex Promise Program is Philips Respironics' commitment that if any patient is noncompliant on any of our Flex technologies, upon return of the unit with Flex technology, we will apply a 100% credit of the corresponding unit purchase price toward the purchase of a Philips Respironics bi-level device with Bi-Flex.

### Bi-level Rescue Program

This program is intended to help patients continue their PAP therapy when they are not able to tolerate standard treatment. Order bi-level monograph (PN 1038429) for more details.

### Extended Warranty Program

The Extended Warranty supports all System One PAP devices and adds three years of coverage to the standard two-year warranty – providing five years of protection on a PAP device or humidifier against unexpected repair costs and shortening the time the device is out of service.

### System One Touch Non-warranty Service Program

For devices not covered under warranty, the System One Touch Non-warranty Service Program offers a simple two-tiered approach to repair costs that helps to accelerate the return of the repaired unit to the patient and reimbursement turnaround times.

## Mask programs

### 30-Day Mask Satisfaction Promise Program

Your patient will be satisfied with our mask or we'll replace it. If a patient discontinues use of a program-approved mask in favor of another mask during the first 30 days of use, we will replace the original mask at no charge. Certain terms and conditions apply.

### Fit for Life Mask Resupply Program

Philips Respironics offers a comprehensive sleep apnea mask solution that helps keep you in touch with your patients through the lifetime of their therapy. By offering a patient interface solution with a mask and resupply service included in the purchase price, Philips Respironics demonstrates a commitment to the long-term success of your patients' therapy and your business. When you choose our masks, your patients will receive access to either our medSage or Encore Resupply services. The Fit for Life solution benefits both your patients and you by increasing patient compliance and satisfaction and helping you develop a standard, quality-driven approach to meeting the payer and industry requirements that are making in-house resupply programs more challenging.



## Ventilation programs

### **Ventilation Promise Program**

Progressive diseases like COPD and ALS have one thing in common – a patient's ventilator needs will change over time. We want to make it easy for you to transition your patients to higher level ventilator devices as their needs become more challenging. To enable you to easily transition patients to other devices as their diseases progress, the Ventilation Promise Program will provide a credit equal to the purchase price of any Philips Respironics' System One device to be applied toward the purchase of a BiPAP AVAPS or Trilogy ventilator. Certain terms and conditions apply.

### **Partners in Ventilation (PIV) Program**

This program has identified the needs of homecare ventilation providers and designed a set of tools with one goal – to simplify your business so you can focus on what matters most – your patients. As a program member you will have access to ventilation-related resources developed by Philips Respironics. And your PIV membership is free.

## Customer programs

### **Partners in Compliance Management**

The Encore Anywhere compliance data management system makes patient data available easily and securely. Additional compliance management resources can be found at: [www.sleepapnea.com](http://www.sleepapnea.com). It's all part of our solution to help you meet your business goals, enhance staff efficiency, and coordinate patient care more quickly and easily.

### **Partners in Training**

The Partners in Training web site is a sleep apnea awareness resource designed to help homecare providers, clinicians, and physicians better understand how to identify, treat and manage sleep apnea, while helping to improve business/practice potential.

### **Cognitive Behavioral Therapy for Insomnia (CBTI)**

This program provides sleep labs with a means to implement a behavioral sleep medicine program or enhance their current offerings without incurring the development cost and time of doing it themselves. Additionally, it enables them to implement an objective tool (Actiwatch) for therapeutic assessment as an adjunct to traditional, subjective diaries.

### **Know How Webinars**

The Philips Respironics Know How program is a series of live and on-demand webinars that help to address key business and clinical challenges and opportunities faced by homecare providers and healthcare professionals. Visit <http://knowhow.respironics.com> to learn more and to register for the various sessions.



# Medical education

For respiratory therapists, sleep technicians and nurses, Philips Respironics CEU-accredited programs target respiratory drug delivery, noninvasive and invasive ventilation, oxygen and sleep therapy and are focused on your continuing education requirements.

To learn about additional training opportunities relating to sleep apnea, visit our Partners in Training web site at [www.pit.respironics.com](http://www.pit.respironics.com)

Not all products are available in all countries. Check with Philips Respironics customer service or your local distributor.



# Sleep therapy presentations



## **The mask experience**

2 CEUs, AARC, AAST

A lecture, followed by a hands-on session to learn how to choose the right mask for the right patient

- Understand the challenges the OSA patient faces in order to be compliant to his or her CPAP therapy
- Understand why your patient needs to be compliant to CPAP therapy
- Understand how to choose the right mask to help your patient be compliant to CPAP therapy
- Understand the types of masks
- Understand the significance of resupply

## **Fundamentals of obstructive sleep apnea**

1 CEU, AARC, AAST, CCMC

Understand the pathology and prevalence of OSA and learn how to recommend appropriate diagnostic and treatment options

- Understand the types of Sleep Disordered Breathing (SDB)
- Understand the prevalence of SDB
- Recommend appropriate diagnosis and screening procedures
- Discuss various treatment and follow-up management options

### Tools for a successful sleep therapy compliance program

1 CEU, AARC, AAST

This session will cover the challenges of running a successful Sleep Therapy Compliance Program and provides caregivers with the clinical tools/techniques to be successful

- Provide overview of PAP compliance
- Understand consequences of non-compliance
- Explore techniques, tools, and technologies to facilitate compliance

### Advanced algorithms for SDB

1.5 CEUs, AARC, AAST

Learn about leak and triggering technology, understand the CPAP breathing algorithms, and identify the breathing algorithms for periodic breathing and respiratory insufficiency

- Review Auto-EPAP, autoSV, AVAPS, Auto-Trak, C-Flex
- Emphasis on Philips Respironics technology

### Titration guidelines for CPAP, bi-level, servo ventilation, and AVAPS

1.5 CEUs, AARC, AAST

List the titration goals for OSA management, understand the suggested protocols for titrations and review the titration process

- Suggest a titration procedure for basic CPAP
- Suggest a titration procedure for autoSV

### Maximizing PAP therapy through advanced techniques: ASV and AVAPS

1 CEU, AARC, AAST

This session will help you identify which of your patients will benefit from advanced therapies

- Explain when to use servo ventilation
- Explain when to use AVAPS

### Data management of the SDB patient

1 CEU, AARC, AAST

Review the management and follow-up of PAP therapy through the use of the Encore data management software and learn how to interpret the reports and waveforms and how this combines with a self-management system

- Describe the management of a patient through the use of data management software
- Recognize that compliance leads to positive health outcomes
- Discuss an overview of compliance technologies in healthcare
- Describe tools and technologies to facilitate PAP compliance

### Home sleep testing

1 CEU, AARC, AAST

A review of the guidelines for home sleep testing and which patient is best suited for diagnosis and titration in the home

- Identify the predisposing factors for OSA
- Describe the co-morbidities that require a full in-lab sleep study
- Identify good candidates for portable monitoring, pre-screening questionnaires
- Describe equipment used to perform portable monitoring



# Respiratory drug delivery presentations



## **Pediatric inhalation devices**

1 CEU, AARC

Presentation will cover pediatric aerosol delivery

- Explain the rationale of use of commonly inhaled drugs that are used in the management of pediatric respiratory disease
- Increase the level of understanding for pediatric aerosol therapy
- Explain the benefits and shortcomings of devices that are used for aerosol therapy for pediatric patients

## **Current aerosol therapy: theory and practice**

1 CEU, AARC, CCMC

Reviews basic respiratory anatomy and physiology and respiratory diseases, provides an overview of the history of inhalation drug delivery devices, and defines the design principles of these devices

- Review of respiratory anatomy and physiology
- Understand basic principles of an aerosol device
- Define and differentiate the following respiratory drug delivery devices: dry powdered inhalers, metered dose inhalers, and jet nebulizers



# Home respiratory care presentations



## **Portable home oxygen equipment matching the patient to the proper device**

**1 CEU, AARC**

This presentation will help participants gain a better understanding of the clinical and operational differences between stationary and portable oxygen equipment and will detail methods for matching portable oxygen equipment to specific patient needs

- Identify methods to match portable oxygen equipment to patient needs
- Discuss the clinical and operational differences among portable oxygen devices/equipment
- Explain how oxygen-conserving device technology applies to portable oxygen equipment

## **Trilogy ventilator workshop**

**2 CEUs, AARC**

Learn more about the operation of the Trilogy100 and Trilogy200 life support ventilators

- Review the circuit options, breath types, modes, alarms, features, and troubleshooting

### CoughAssist T70 approach to secretion management

1 CEU,AARC

This workshop educates respiratory therapists and nurses on the operation of the CoughAssist T70 device and will cover indications, contraindications, operation procedure recommended treatment practice and review of a clinical case study

- Understand the indications, contraindications and warnings/cautions for the use of CoughAssist T70
- Identify the advantages of the noninvasive use of CoughAssist T70 over invasive airway clearance
- Review a clinical case study
- Learn the controls, connectors, and visual indicators of the device
- Discuss the recommended treatment plan
- Understand cleaning and maintenance of the CoughAssist T70 device

### The continuum of products and solutions for the COPD patient

1 CEU,AARC

Participants will receive a brief review of COPD staging as we discuss how Philips Respironics' comprehensive product portfolio can help you to manage your COPD patients throughout the disease progression

- Review COPD
- Discuss current practices for treatment and components of care
- Review COPD staging
- Identify the Philips Respironics products and solutions as they relate to the progression of COPD
- Present tools which will help your DME to grow and thrive in today's market



### Respiratory advances of mouth piece ventilation

1 CEU,AARC

This course will cover an additional option to provide noninvasive ventilation for the chronic respiratory insufficiency patient

- Identify options for 24-hour ventilation without the need for tracheostomy
- Discuss how to prevent the pulmonary complications related to tracheostomy
- Describe infection reduction–granuloma formation chronic bacterial colonization
- Explain cough generation from a closed airway

Philips Respironics-sponsored online courses can be found at [www.ganESCO.com](http://www.ganESCO.com).



**Philips Healthcare is part of  
Royal Philips**

**How to reach us**

[www.philips.com/healthcare](http://www.philips.com/healthcare)  
[healthcare@philips.com](mailto:healthcare@philips.com)

Asia  
+49 7031 463 2254

Europe, Middle East, Africa  
+49 7031 463 2254

Latin America  
+55 11 2125 0744

North America  
+1 425 487 7000  
800 285 5585 (toll free, US only)

Philips Respironics  
1010 Murry Ridge Lane  
Murrysville, PA 15668

Customer Service  
+1 724 387 4000  
800 345 6443 (toll free, US only)

Philips Respironics Asia Pacific  
+65 6882 5282

Philips Respironics Australia  
+61 (2) 9947 0440  
1300 766 488 (toll free Australia only)

Philips Respironics Benelux  
+33 1 47 28 30 83

Philips Respironics China  
+86 400 828 6665  
+86 800 828 6665

Philips Respironics Deutschland  
+49 8152 93 06 0

Philips Respironics France  
+33 2 51 89 36 00

Philips Respironics Iberia  
+34 91 566 90 00

Philips Respironics Italy  
+39 039 203 1

Philips Respironics Sweden  
+46 8 120 45 900

Philips Respironics Switzerland  
+41 6 27 45 17 50

Philips Respironics United Kingdom  
+44 800 1300 845

[www.philips.com/respironics](http://www.philips.com/respironics)

Please visit [www.philips.com/respironics](http://www.philips.com/respironics)



©2014 Koninklijke Philips N.V.  
All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

CAUTION: US federal law restricts these devices to sale by or on the order of a physician.

edoc JH 06/06/14 MCI 4106218